


evercare[®] Ophthalmic set, incision film

REF: 1280-01

evercare[®] Ophthalmic set with integrated incision film, sterile contains an impermeable eye drape for efficient infection control. It has an integrated pouch for fluid handling and a flexible wire to maintain wanted shape in fluid collection pouch. Also equipped with incise film in the aperture field to provide full access.

Key Features

- ▶ Incision film
- ▶ Fluid collection pouch
- ▶ Main drape complies to EN 13795-1 High performance critical product area requirement




REF	DESCRIPTION	SIZE	Pack Count dispenser/carton
1280-01	Ophthalmic set with integrated incision film, sterile	See Material table	7 / 28


Sales Unit: Transport Carton

Regulations

The product complies with legal requirements for medical devices. EO sterilized according to valid version of ISO 11135 standard.




2460




STERILE


EO


MD




Not made with natural rubber latex












15 °C - 25 °C



2

STERILIZATION





REACH regulation:	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
Biological evaluation:	EN ISO 10993-01 EN ISO 10993-5 EN ISO 10993-7 EN ISO 10993-10
Compliance with legal requirements:	Medical Device Regulation (EU) 2017/745, MDR class Is
Compliance with product standards:	EN13795-1
Quality standards:	EN ISO 13485
Label information and packaging compliance:	EN ISO 15223-1, EN ISO 20417, EN ISO 11607-1

Material

Product Part	Raw Material & Characteristics
1 Hand towel	30 x 34 cm, Pulp, natural white, thread reinforced
1 Ophthalmic drape, with incision film	150 x 140 cm, laminate: PP nonwoven / PE-film blue Fluid collection pouch: PE/EPA-film, transparent Metal wire coated with PE, blue Incision film: adhesive coated polyethylene film
1 Instrument table cover	150 x 140 cm, laminate: PP nonwoven / PE-film, blue

Page 2

2023-12-05

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur. Do not use if package is damaged.
Storage Recommendation	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 years

Packing Information

Package	Material	Size (L x W x H) mm
Single Pack	Medical grade paper, PE/PET film	250 x 390
Dispenser Pack	TCF or ECF quality	285 x 242 x 287
Transport Carton	TCF or ECF quality	500 x 310 x 590
Packing Tape	Transparent	N/A

Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
1280-01	06438129110333	06438129310337	06438129510331

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.