

evercare[®] Absorbent waterproof adhesive dressing, sterile

REF: 4655, 4657, 4660, 4665, 4670, 4675, 4680, 4685

evercare[®] Absorbent waterproof adhesive dressing, sterile is used for surgical wound treatment and exudate absorption. The absorbent pad is finished by semi-permeable net, allowing to absorb significant quantity of exudate coming from the wound while performing non-adherent properties at the same time. High retention rate allows to keep exudate inside absorbent pad even under external pressure. The dressing is air- and moisture vapor permeable. It can easily be applied on conical and round part of the patient's body. The PU foil is semipermeable and protects against moisture and contamination.

Key Features

- ▶ Soft and flexible
- ▶ High breathability
- ▶ For low to moderately exuding wounds
- ▶ Acrylic based adhesive
- ▶ Easy-to-use unit pouch
- ▶ Non-woven material covered by PU foil
- ▶ Showerable




REF	DESCRIPTION	SIZE	Pack Count dispenser/carton
4655	Absorbent waterproof adhesive dressing with pad, sterile	5x7 cm; pad 2,5x4 cm	50 / 2400
4657	Absorbent waterproof adhesive dressing with pad, sterile	10x6 cm; pad 5x3 cm	50 / 1600
4660	Absorbent waterproof adhesive dressing with pad, sterile	10x10 cm; pad 4,5x4,5 cm	50 / 1200
4665	Absorbent waterproof adhesive dressing with pad, sterile	10x15 cm; pad 5x11 cm	50 / 600
4670	Absorbent waterproof adhesive dressing with pad, sterile	10x20 cm; pad 5x15,5 cm	30 / 360
4675	Absorbent waterproof adhesive dressing with pad, sterile	10x25 cm; pad 5x20 cm	30 / 360
4680	Absorbent waterproof adhesive dressing with pad, sterile	10x30 cm; pad 5x25 cm	30 / 360
4685	Absorbent waterproof adhesive dressing with pad, sterile	10x35 cm; pad 5x30 cm	30 / 360

Sales Unit: Transport Carton

Regulations

The product complies with legal requirements for medical devices. EO sterilized according to valid version of ISO 11135 standard.



REACH regulation: The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).

Biological evaluation: EN ISO 10993-5
EN ISO 10993-7
EN ISO 10993-10
EN ISO 10993-11

Compliance with legal requirements: Medical Device Regulation (EU) 2017/745, MDR class IIa

Compliance with product standards: EN ISO 11135:2014/A1:2019, EN ISO 11737-1, EN 13726-2

Quality standards: EN ISO 13485

Label information and packaging compliance: EN ISO 15223-1, EN ISO 20417, EN ISO 11607-1

Material

Product Part	Raw Material & Characteristics
Carrier	PU breathable and waterproof + 100% polyester nonwoven, water repellent, white
Coating	Pressure sensitive glue, thermosetting acrylic solution polymer, latex free
Absorbent pad	70 % viscose, 30 % polyester and non-adherent porous PE film, white
Release paper	Coated paper, white

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space between 15 - 25 °C. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	3 years

Packing Information

Package	Material	Size (L x W x H) mm
Single Pack	Medical grade paper, PP/PE film	4655: 72 x 113
		4657: 135 x 124
		4660: 135 x 160
		4665: 135 x 220
		4670: 135 x 262
		4675: 135 x 313
		4680: 135 x 361
4685: 135 x 413		
Dispenser Pack	Carton, TCF or ECF quality	4655: 95 x 74 x 115
		4657: 140 x 95 x 130
		4660: 140 x 95 x 165
		4665: 140 x 95 x 225
		4670: 140 x 65 x 265
		4675: 140 x 65 x 316
		4680: 140 x 65 x 365
4685: 140 x 65 x 415		
Transport Carton	Carton, TCF or ECF quality	4655: 320 x 600 x 260
		4657: 400 x 600 x 280
		4660: 300 x 600 x 360
		4665: 600 x 300 x 360
		4670: 285 x 440 x 285
		4675: 285 x 440 x 335
		4680: 285 x 440 x 385
4685: 285 x 440 x 415		
Packing Tape	PVC-free, Transparent	N/A

Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
4655	06438129115383	06438129315387	06438129515381
4657	06438129118568	06438129318562	06438129518566
4660	06438129115390	06438129315394	06438129515398
4665	06438129115406	06438129315400	06438129515404
4670	06438129115413	06438129315417	06438129515411
4675	06438129115420	06438129315424	06438129515428
4680	06438129115437	06438129315431	06438129515435
4685	06438129115444	06438129315448	06438129515442

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.