

evercare[®] MediKit Dressing Pack, sterile

REF: 53693

evercare[®] Dressing Pack, sterile is a convenient set that support efficiency and can help reduce preparation time, including components needed for specific interventions.

Key Features

- ▶ 1 Plastic tray, 3 compartments, 150x76x18 mm
- ▶ 8 Non woven sponges diam. 30 mm
- ▶ 1 Forceps, 14 cm
- ▶ 1 Laminated paper, 40x38 cm




REF	DESCRIPTION	SIZE	Pack Count dispenser/carton
53693	Dressing Pack, sterile	N/A	25 / 200

Sales Unit: Transport Carton

Regulations

The product complies with legal requirements for medical devices. EO sterilized according to valid version of ISO 11135 standard.



REACH regulation:	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
Biological evaluation:	EN ISO 10993-5 EN ISO 10993-7 EN ISO 10993-10
Compliance with legal requirements:	Medical Devices Directive 93/42/EEC, MDD class Is
Compliance with product standards:	EN ISO 11135
Quality standards:	EN ISO 13485
Label information and packaging compliance:	EN ISO 15223-1, EN ISO 20417

Material

Product Part	Raw Material & Characteristics
2 Plastic tray, 3 compartments, 150x76x18 mm	PS
8 Non woven sponges diam. 30 mm	Rayon/PES, white
1 Forceps, 14 cm	PS
1 Laminated paper, 40x38 cm	Paper/PE, blue/white

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space between 0 - 25 °C. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 years

Packing Information

Package	Material	Size (L x W x H) mm
Single Pack	Medical grade paper, Plastic film	133 x 222
Dispenser Pack	Carton, TCF or ECF quality	280 x 180 x 200
Transport Carton	Carton, TCF or ECF quality	590 x 390 x 430
Packing Tape	PVC-free, Transparent	N/A

Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
53693	06937935310093	06937935340090	06937935370097

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.