

## evercare® Elastic fixation bandage Crepe

REF: 4161-01, 4162-01, 4163-01, 4164-01, 4165-01

evercare® Elastic fixation bandage Crepe is stretchable with soft and porous structure. Due to its air permeability and softness the bandages are easy to apply, conformable and comfortable to use. Thanks to woven-edge weaving technique and careful selection of raw material, the bandages may be used for fixation of most dressings, especially for joints and round body parts, and small medical devices to the patient's body. Package designed for easy product identification.

### Key Features

- ▶ Fully woven edges to minimize loose threads
- ▶ Fixation of all types of dressings, especially for joints and round body parts
- ▶ Soft and comfortable
- ▶ Crepe structure for extra non-slip function
- ▶ 90% stretchability
- ▶ 100% polyester




REF	DESCRIPTION	SIZE	Pack Count dispenser/carton
4161-01	Elastic fixation bandage Crepe, polyester, White	4 cm x 4 m	20 / 600
4162-01	Elastic fixation bandage Crepe, polyester, White	6 cm x 4 m	20 / 560
4163-01	Elastic fixation bandage Crepe, polyester, White	8 cm x 4 m	20 / 480
4164-01	Elastic fixation bandage Crepe, polyester, White	10 cm x 4 m	20 / 400
4165-01	Elastic fixation bandage Crepe, polyester, White	12 cm x 4 m	20 / 320

Sales Unit: Transport Carton

## Regulations

The product complies with legal requirements for medical devices.



<b>REACH regulation:</b>	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
<b>Biological evaluation:</b>	EN ISO 10993-5 EN ISO 10993-10
<b>Compliance with legal requirements:</b>	Medical Device Regulation (EU) 2017/745, MDR class I
<b>Compliance with product standards:</b>	DIN 61632, DIN 61634, EN 1049-2
<b>Quality standards:</b>	EN ISO 13485
<b>Label information and packaging compliance:</b>	EN ISO 15223-1, EN ISO 20417

## Material

Product Part	Raw Material & Characteristics
<b>Bandage</b>	100% polyester, white, Unstretched length 2 m, Stretched length 4 m, Elongation approx. 90%, Elasticity approx. 80%

## Warnings & Recommendations

<b>Warnings</b>	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
<b>Storage Recommendation</b>	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
<b>Disposal Recommendation</b>	Dispose in accordance with local regulations.
<b>Shelf Life</b>	5 years

## Packing Information

Package	Material	Size ( L x W x H ) mm
<b>Single Pack</b>	N/A	N/A
<b>Dispenser Pack</b>	Carton, TCF or ECF quality	4161-01: 165 x 45 x 135 4162-01: 165 x 65 x 135 4163-01: 165 x 85 x 135 4164-01: 165 x 105 x 135 4165-01: 165 x 125 x 135
<b>Transport Carton</b>	Carton, TCF or ECF quality	4161-01: 425 x 350 x 245 4162-01: 485 x 350 x 290 4163-01: 540 x 350 x 290 4164-01: 550 x 350 x 290 4165-01: 520 x 350 x 290
<b>Packing Tape</b>	PVC-free, Transparent	N/A

## Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
<b>4161-01</b>	N/A	06438129320992	06438129520996
<b>4162-01</b>	N/A	06438129321005	06438129521009
<b>4163-01</b>	N/A	06438129321012	06438129521016
<b>4164-01</b>	N/A	06438129321029	06438129521023
<b>4165-01</b>	N/A	06438129321036	06438129521030

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.