



SELEFA® Triangular gauze sponges, sterile

REF: 222241

SELEFA® Triangular gauze sponges, sterile are highly absorbent sponges for a wide area of use such as: skin cleansing, wound cleansing, wiping, fluid and moisture absorption, protection, disinfection etc. Good absorbent properties and in compliance with EN 14079 to assure high quality cotton gauze. The compact folding keeps the product shape stable. Gauze sponges in these small sizes are suitable for various applications where sterility is also required, as cleaning of small wounds and abrasions, minor operation preppings, debris removal etc.

Key Features

- ▶ 100% cotton gauze
- ▶ Environmental friendly bleached (no chlorine)
- ▶ Folded edges prevent loose threads
- ▶ Sterilized using EO
- ▶ 5-pack
- ▶ Easy-to-use unit pouch
- ▶ 17 threads per cm²




REF	DESCRIPTION	SIZE	Pack Count dispenser/carton
222241	Triangular gauze sponges, 5-pack, sterile, White	30x30x30 mm	250 / 5000

Sales Unit: Transport Carton

Regulations

The product complies with legal requirements for medical devices. EO sterilized according to valid version of ISO 11135 standard.



REACH regulation:	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
Biological evaluation:	EN ISO 10993-5 EN ISO 10993-7 EN ISO 10993-10
Compliance with legal requirements:	Medical Device Regulation (EU) 2017/745, MDR class Is
Compliance with product standards:	EN 14079, EN ISO 11135
Quality standards:	EN ISO 13485
Label information and packaging compliance:	EN ISO 15223-1, EN ISO 20417, EN ISO 11607-1

Material

Product Part	Raw Material & Characteristics
Main material	100 % cotton gauze, white

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space between 15 - 25 °C. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 years

Packing Information

Package	Material	Size (L x W x H) mm
Single Pack	Medical grade paper, PP/PE film	
Dispenser Pack	Carton, TCF or ECF quality	237 x 100 x 116
Transport Carton	Carton, TCF or ECF quality	595 x 245 x 415
Packing Tape	PVC-free, Transparent	N/A

Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
222241	06438129154108	06438129354102	06438129554106

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.