

SELEFA® Gauze swabs

REF: 22219512, 22219513, 22219523, 22219524, 22219534, 22219535, 22219545, 22219546

SELEFA® Gauze swabs are primarily intended to cover or clean and absorb blood and exudate from wounds and abrasions. Can also be used for non-invasive medical procedures, such as wound padding, drying, protection, skin and wound cleansing or as secondary dressing. Good absorbent properties and in compliance with EN 14079 to assure high quality cotton gauze.

Key Features

- ▶ 100% cotton gauze
- ▶ Environmental friendly bleached (no chlorine)
- ▶ 17 threads per cm²



REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
22219512	Gauze swabs, non-sterile, White	5x5 cm 12-ply	100 / 5000
22219513	Gauze swabs, non-sterile, White	5x5 cm 8-ply	100 / 5000
22219523	Gauze swabs, non-sterile, White	7.5x7.5 cm 12-ply	100 / 5000
22219524	Gauze swabs, non-sterile, White	7.5x7.5 cm 8-ply	100 / 5000
22219534	Gauze swabs, non-sterile, White	10x10 cm 12-ply	100 / 5000
22219535	Gauze swabs, non-sterile, White	10x10 cm 8-ply	100 / 5000
22219545	Gauze swabs, non-sterile, White	10x20 cm 12-ply	100 / 2000
22219546	Gauze swabs, non-sterile, White	10x20 cm 8-ply	100 / 3000

Sales Unit: Transport Carton

Regulations

The product complies with legal requirements for medical devices.



REACH regulation:	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
Biological evaluation:	EN ISO 10993-5 EN ISO 10993-10
Compliance with legal requirements:	Medical Device Regulation (EU) 2017/745, MDR Class I
Compliance with product standards:	EN 14079
Quality standards:	EN ISO 13485
Label information and packaging compliance:	EN ISO 15223-1, EN ISO 20417

Material

Product Part	Raw Material & Characteristics
Main material	100 % cotton gauze, white

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 years

Packaging Information

Package	Material	Size (L x W x H) mm
Single Pack	N/A	N/A
Dispenser Pack	Medical grade paper	22219512: 50 x 50 x 170 22219513: 200 x 100 x 170 22219523: 50 x 50 x 160 22219524: 50 x 50 x 120 22219534: 75 x 75 x 170 22219535: 75 x 75 x 120 22219545: 100 x 100 x 170 22219546: 100 x 100 x 130
Transport Carton	Carton	22219512: 260 x 260 x 350 22219513: 415 x 520 x 355 22219523: 265 x 270 x 335 22219524: 265 x 270 x 255 22219534: 390 x 395 x 355 22219535: 390 x 395 x 255 22219545: 515 x 520 x 355 22219546: 515 x 520 x 275
Packing Tape	Transparent	N/A

Barcodes / UDI-DI

REF	Single Pack	Dispenser Pack	Transport Carton
22219512	N/A	06438129351835	06438129551839
22219513	N/A	06438129352252	06438129552256
22219523	N/A	06438129351842	06438129551846
22219524	N/A	06438129352269	06438129552263
22219534	N/A	06438129351859	06438129551853
22219535	N/A	06438129352276	06438129552270
22219545	N/A	06438129351866	06438129551860
22219546	N/A	06438129352283	06438129552287

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.