

evercare[®] inLine 3-way stopcock, Microtubing, Extension line, PVC-free

REF: 10699-01

evercare[®] inLine 3-way stopcocks are a and broad range of product for different needs. 3-way stopcock with extension tubing in different legths enabling connections to two sources of infusion and/or injection or a combination of both. Stopcock are sterile, single use products that are used to conduct intravenous fluids as an extension of the infusion or perfusion device, by enabling connections to two sources of infusion and an infusion device with vascular access.

Key Features

- ▶ Lipid resistant
- ▶ 45° click feature
- ▶ Products can withstand upto 2 bar liquid pressure
- ▶ Hydrophobic filter prevents leakage while priming
- ▶ Rotating luer lock




REF	DESCRIPTION	SIZE - Tube Inner/outer, L(length)	PRIMING VOLUME(ml) Approx	Pack Count dispenser/carton
10699-01	3-way stopcock, Microtubing, Extension line, PVC-free	Ø 1/3.1 mm, L 15cm	0.4	25 / 100

Sales Unit: Transport Carton

Regulations

The product complies with legal requirements for medical devices. EO sterilized according to valid version of ISO 11135 standard.



REACH regulation:

The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).

Biological evaluation:

EN ISO 10993-01

Compliance with legal requirements:

Medical Device Regulation (EU) 2017/745, MDR Class IIa

Compliance with product standards:

EN 8536-9, EN ISO 80369-7

Quality standards:

EN ISO 13485

Label information and packaging compliance:

EN ISO 15223-1, EN ISO 20417

Material

Product Part	Raw Material & Characteristics
Male luer cap	HDPE, white
3-way stopcock	Transparent channel PC, lipid resistant, white handle HDPE
Tube	Polyolefin alloy, PVC-free, soft, clear, kink resistant
Rotating luer lock	ABS & PC, transparent
Purge filter cap	HDPE, light blue

Warnings & Recommendations

Warnings	<p>The use of this product is restricted to a qualified doctor or nurse.</p> <p>Read the instructions before use. The device should be used in accordance with the instructions.</p> <p>The device is guaranteed sterile and pyrogen-free if the individual package has not been opened or damaged. Do not use if package is damaged or if protective caps are loose or missing.</p> <p>Do not over-tighten a luer beyond the friction fit as this may damage both the luer and another connecting device.</p> <p>Do not use any instruments for clamping the tubing or luer tightening.</p> <p>Re-use may lead to infection or other illness/injury.</p> <p>Do not leave open ports/hubs exposed. Replace sterile end caps, as necessary.</p>
Storage Recommendation	Recommended storage in clean, dry space between 15 - 25 °C. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 years

Packing Information

Package	Material	Size (L x W x H) mm
Single Pack	Medical grade paper, PP/PE film	235 x 100
Dispenser Pack	PP	305 x 460 x 0
Transport Carton	TCF quality, Carton	400 x 295 x 287
Packing Tape	PVC-free, Transparent	N/A

Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
10699-01	06438129121575	N/A	06438129521573

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.