




evercare<sup>®</sup> MediKit Round bowl 120 ml, 5 gauze swabs, sterile  
REF: 191020


evercare<sup>®</sup> Round bowl 120 ml, 5 gauze swabs, sterile is a convenient set that support efficiency and can help reduce preparation time, including components needed for specific interventions.

Key Features

- ▶ 1 Cup, 120 ml
- ▶ 5 Gauze swabs, 5x5 cm



5 Gauze swabs, 5 × 5 cm



1 Cup, 120 ml

REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
191020	Round bowl 120 ml, 5 gauze swabs, sterile	N/A	105 / 210

Sales Unit: Transport Carton

## Regulations

The product complies with legal requirements for medical devices. EO sterilized according to valid version of ISO 11135 standard.



REACH regulation:	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
Biological evaluation:	EN ISO 10993-5 EN ISO 10993-7 EN ISO 10993-10
Compliance with legal requirements:	Medical Devices Directive 93/42/EEC, MDD Class Is
Compliance with product standards:	EN 14079, EN ISO 11135
Quality standards:	EN ISO 13485
Label information and packaging compliance:	EN ISO 15223-1, EN ISO 20417, EN ISO 11607-1

## Material

Product Part	Raw Material & Characteristics
1 Cup, 120 ml	PP, transparent
5 Gauze swabs, 5x5 cm	Cotton, 13 threads/cm², white

## Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space between 0 - 25 °C. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 years

## Packaging Information

Package	Material	Size ( L x W x H ) mm
Single Pack	Medical grade paper, Plastic film	131 x 156
Dispenser Pack	Carton, TCF or ECF quality	378 x 283 x 410
Transport Carton	Carton, TCF or ECF quality	590 x 390 x 430
Packing Tape	PVC-free, Transparent	N/A

## Barcodes / UDI-DI

REF	Single Pack	Dispenser Pack	Transport Carton
191020	06937935320207	06937935330169	06937935360203

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.