

Uterine Aspiration Tubing

Product Data Sheet

Description	A range of flexible transparent tubing offering a convenient and adaptable system for uterine aspiration.
Intended Use	A sterile, semi-rigid or rigid, plastic tube designed for insertion into the uterus. It is frequently connected directly or through an external tube or catheter to an irrigator and/or aspirator for irrigation, infusion, or withdrawal of substances. It may also provide a port through which instruments and/or probes are passed to perform a variety of gynaecologic/obstetric diagnostic, therapeutic, and surgical procedures. This is a single-use device. (GMDN 32655)
Duration of Use	Normally intended for continuous use of less than 60 minutes, but will be patient specific depending on duration of the surgical procedure.
Materials	PVC and K-Resin. All materials used are DEHP and LATEX FREE.
Product Information	Special wide bore, thick walled, high vacuum PVC tubing, with a Flexi-cuff to provide extra flexibility at the crucial point. Attached to rotary swivel handle featuring easy vacuum control. Flexi-Cuff provides extra flexibility enabling ease-of-use and potentially enhancing techniques. Leak proof seals; Easily connects to curettes; Male connector in each pack for connection to suction machines where required; Sterile, double wrapped, individual peel pouches.
Quality System	ISO 13485
EU Classification	Class IIa
Sterilised by	Ethylene Oxide;
Shelf Life	5 years
Storage & Safety	Toxicity: Manufacture from biocompatible materials Storage: Store in a clean, cool, dry place away from direct sunlight



Included Male Connector:



Tubing Length:	2m or 3m
Tubing Diameter:	11.5mm (ID) 16mm (OD)
Tubing Shore Hardness:	75A
Flexi-Cuff Shore Hardness:	62A

Product Code	Description	Tubing Length	Qty per Box
AEKHS3000	Wide Bore Tubing + Rotary Handle, loose male connector	2.0m	10
AEKHS3000M	Wide Bore Tubing + Rotary Handle, attached male connector	2.0m	10
AEKHS3030	Wide Bore Tubing + Rotary Handle, loose male connector	3.0m	10



 GCMEDICA ENTERPRISE LTD., (WUXI)
 Loujin Industrial Park, Shuofang, Wuxi 214143, Jiangsu, China
 Luxus Lebenswelt GmbH
 Kochstr.1, 47877, Willich, Germany

© 2021 Pennine Healthcare or one of its subsidiaries. All rights reserved.

The information given in this sheet is not intended to be exhaustive.

The information contained in this sheet is liable to modification from time to time in the light of experience and our policy of continuous product development.