

evercare[®] Adhesive Plasters nonwoven

REF: CMT352177, CMT352184, CMT352191, CMT352207

evercare[®] Adhesive plaster nonwoven is a dressing intended to protect smaller wounds from impurities in outer surrounding. The absorbent patch attached to the plaster absorbs exudate from the covered wound. The plaster material is soft and comfortable for the skin. The plaster is typically used for post-injection site protection.

Key Features

- ▶ Soft, high breathable nonwoven plaster
- ▶ Solvent-free adhesive
- ▶ Non-sticking absorbent pad
- ▶ Two different shapes: stripe and rounded
- ▶ Individual pouch for each plaster




REF	DESCRIPTION	SIZE	Pack Count dispenser/carton
CMT352177	Adhesive Plaster, nonwoven, beige, stripe; Producer Item No. CMT352177, Beige	2,4x7cm	200 / 10200
CMT352184	Adhesive Plaster, nonwoven, beige, stripe; Producer Item No. CMT352184, White	1,2x4cm	200 / 10200
CMT352191	Adhesive Plaster, nonwoven, beige, stripe; Producer Item No. CMT352191, Beige	3,8x3,8cm	200 / 10200
CMT352207	Adhesive Plaster, nonwoven, beige, stripe; Producer Item No. CMT352207, Beige	2x4cm	200 / 10000

Sales Unit: Transport Carton

Regulations

The product complies with legal requirements for medical devices.



REACH regulation:	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
Biological evaluation:	EN ISO 10993-01 EN ISO 10993-5 EN ISO 10993-10
Compliance with legal requirements:	Medical Device Regulation (EU) 2017/745, MDR class I
Compliance with product standards:	EN 14079, ASTM D3330
Quality standards:	EN ISO 13485
Label information and packaging compliance:	EN ISO 15223-1, EN ISO 20417

Material

Product Part	Raw Material & Characteristics
Main material	Nonwoven, 100 Polyester, beige
Adhesive	Hot-melt acrylic glue, UV cross linked
Absorbent pad	Viscose/PP/PE, layered with non-stick LDPE
Peel-paper	Silicone paper, one sided

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 years

Packing Information

Package	Material	Size (L x W x H) mm
Single Pack	N/A	N/A
Dispenser Pack	Corrugated cardboard, TCF or ECF quality	93 x 165 x 35
Transport Carton	Corrugated cardboard, TCF or ECF quality	580 x 185 x 300
Packing Tape	Transparent, PVC-free	N/A

Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
CMT352177	N/A	N/A	N/A
CMT352184	N/A	N/A	N/A
CMT352191	N/A	N/A	N/A
CMT352207	N/A	N/A	N/A

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.