

492**ORD**
OPERATION ROOM DRESSINGS
Xray

Operation dressings

Pensos cirúrgicos

Pansements operatoires

Op-Abdecktuch

Operatieverbanden

Medicazioni per interventi

Apósitos quirúrgicos

Operations textilier



Sterile gauze dressings

- 100% cotton.
- 17 threads/sq.cm (approx. 23g/sq.m.), 20 threads/sq.cm (approx. 27/sq.m.) according to European Pharmacopea and EN 14079.
- X-Ray detectable thread with approx. 60% barium sulphat contents.
- Products are counted and weight controlled to assure maximum efficiency of counting the exact number of products in the package, and packed in the inner package where a double control card is fixed. The inner pack is then packed in an outer peel-open pack and sterilized.
- Inner pack with double control tag + outer peel-open pack + shelf box + transport carton.

Sterile - Double pack - Control tag

Code	Description	Size	Plies	X-Ray	Threads	Pcs./Pack	Pcs./Sh.box
492-006	Gauze swabs	10 x 10cm	12	X	17	5	100
492-007	Gauze swabs	10 x 60cm	4	X	17	5	120
492-008	Gauze swabs	10 x 60cm	8	X	17	5	60
492-010	Laparotomy sponges	30 x 45cm	4	X	20	5	50
492-011	Laparotomy sponges	45 x 70cm	4	X	20	1	20
492-012	Laparotomy sponges	45 x 70cm	4	X	20	5	25
492-013	Laparotomy sponges	10 x 60cm	4	X	17	5	100
492-014	Laparotomy sponges	10 x 60cm	8	X	17	5	60
492-015	Laparotomy sponges	10 x 60cm	12	X	17	5	45