






# evercare<sup>®</sup> MediKit Dressing set, gauze, sterile

REF: 1963-02

evercare<sup>®</sup> Dressing set, gauze, sterile is a compact set containing supplies for dressing change.

## Key Features

- ▶ 5 Round sponges, Ø 30 mm
- ▶ 2 Gauze swabs, 8-ply, 5x5 cm
- ▶ 1 Forceps, 127 mm
- ▶ 1 Tray, 4-compartment , 150 x 76 x 18 mm
- ▶ 1 Utility drape (Tray wrap) 40x38 cm




REF	DESCRIPTION	SIZE	Pack Count dispenser/carton
1963-02	Dressing set, gauze, sterile	N/A	25 / 200

Sales Unit: Transport Carton

## Regulations

The product complies with legal requirements for medical devices. EO sterilized according to valid version of ISO 11135 standard.



**REACH regulation:** The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).

**Biological evaluation:** EN ISO 10993-5  
EN ISO 10993-7  
EN ISO 10993-10

**Compliance with legal requirements:** Medical Devices Directive 93/42/EEC, MDD class Is

**Compliance with product standards:** EN ISO 11135

**Quality standards:** EN ISO 13485

**Label information and packaging compliance:** EN ISO 15223-1, EN ISO 20417, EN ISO 11607-1

## Material

Product Part	Raw Material & Characteristics
5 Round sponges, Ø 30 mm	100 % cotton gauze, 17-thread, elastic ring SEBS/PP
2 Gauze swabs, 8-ply, 5x5 cm	100 % cotton gauze, 17-thread
1 Forceps, 127 mm	Polystyrene, green
1 Tray, 4-compartment , 150 x 76 x 18 mm	Polystyrene
1 Utility drape (Tray wrap) 40x38 cm	Laminate: tissue/pe, green

## Warnings & Recommendations

<b>Warnings</b>	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
<b>Storage Recommendation</b>	Recommended storage in clean, dry space between 0 - 25 °C. Products should be protected from direct sunlight, other intensive light sources and ozone.
<b>Disposal Recommendation</b>	Dispose in accordance with local regulations.
<b>Shelf Life</b>	5 years

## Packing Information

Package	Material	Size ( L x W x H ) mm
<b>Single Pack</b>	Medical grade paper, Plastic film	222 x 133
<b>Dispenser Pack</b>	Carton, TCF or ECF quality	280 x 200 x 180
<b>Transport Carton</b>	Carton, TCF or ECF quality	590 x 430 x 390
<b>Packing Tape</b>	PVC-free, Transparent	N/A

## Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
<b>1963-02</b>	06937935320320	06937935330268	06937935360326

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.