




evercare[®] MediKit Washing set, nonwoven, sterile

REF: 1985-02

evercare[®] Washing set, nonwoven, sterile is a compact sterile set containing round nonwoven sponges for skin disinfection.

Key Features

- ▶ 10 Round nonwoven sponges, diam. 50 mm, packed in double pack



REF	DESCRIPTION	SIZE	Pack Count dispenser/carton
1985-02	Washing set, nonwoven, sterile	N/A	20 / 160

Sales Unit: Transport Carton

Regulations

The product complies with legal requirements for medical devices. EO sterilized according to valid version of ISO 11135 standard.

CE 2460

CE

STERILE EO

MD

LATEX

Not made with natural rubber latex

2

STERILIZE

0 °C

25 °C

REACH regulation:	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
Biological evaluation:	EN ISO 10993-5 EN ISO 10993-7 EN ISO 10993-10
Compliance with legal requirements:	Medical Devices Directive 93/42/EEC, MDD class Is
Compliance with product standards:	EN ISO 11135
Quality standards:	EN ISO 13485
Label information and packaging compliance:	EN ISO 15223-1, EN ISO 20417, EN ISO 11607-1

Material

Product Part	Raw Material & Characteristics
10 Round nonwoven sponges, diam. 50 mm, packed in double pack	Nonwoven, rayon/PES, elastic ring SEBS/PP

Page 2

2023-12-05

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space between 0 - 25 °C. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 years

Packing Information

Package	Material	Size (L x W x H) mm
Single Pack	Medical grade paper, Plastic film	179 x 198
Dispenser Pack	Carton, TCF or ECF quality	280 x 200 x 180
Transport Carton	Carton, TCF or ECF quality	590 x 430 x 390
Packing Tape	PVC-free, Transparent	N/A

Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
1985-02	06937935320405	06937935330343	06937935360401

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.