

evercare[®] Turtleneck poncho

REF: 5026980

evercare[®] Turtleneck poncho with a white knitted neckline for extra warmth and comfort which can be used by both staff and patients when needed to stay warm. One size fits all.

Key Features

- ▶ evercare[®] poncho with turtleneck is comfortable and warm to wear.
- ▶ Can be used for both staff and patients when needed to stay warm.
- ▶ Soft and breathable non-woven material with a very low linting level.
- ▶ The neck cuff has a really good elasticity which allows it to be used by all irrespective of size.



| REF | DESCRIPTION | SIZE | Pack Count dispenser/carton |
|---------|---|------------|--------------------------------|
| 5026980 | Turtleneck poncho with knitted neckline, Blue | 56 x 46 cm | 40 / 400 |

Sales Unit: Transport Carton

Regulations

N/A



Not made with
natural rubber latex



REACH regulation:

The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).

Material

| Product Part | Raw Material & Characteristics |
|------------------|---|
| Nonwoven | Nonwoven PP SMS, blue |
| Knitted neckline | Knitted neckline PES rib knitted, white |
| Sewing thread | Sewing thread Elastic Polyester thread, white |

Warnings & Recommendations

| Warnings | |
|-------------------------|---|
| Storage Recommendation | Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone. |
| Disposal Recommendation | Dispose in accordance with local regulations. |
| Shelf Life | 5 years |

Packing Information

| Package | Material | Size (L x W x H) mm |
|------------------|--------------------|-----------------------|
| Single Pack | N/A | N/A |
| Dispenser Pack | TCF or ECF quality | 250 x 170 x 270 |
| Transport Carton | TCF or ECF quality | 560 x 360 x 520 |
| Packing Tape | Transparent | N/A |

Barcodes / UDI-DI

| REF | Single Pack | Dispenser | Transport Carton |
|---------|-------------|---------------|------------------|
| 5026980 | N/A | 6438129325188 | 6438129525182 |

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.