



evercare SUPRA Standard OP-Gowns, sterile

REF: 259011, 259021, 259031, 259041, 259051, 259061, 259071, 259081, 259091, 259101, 259111, 259121, 259131, 259141

evercare® SUPRA OP-Gown is a soft SMS gown designed to provide ultimate protection and super comfort. The gown is suitable for procedures with low to medium to high quantities of fluids and high comfort demands. SUPRA OP-Gown material has antistatic and alcohol repellent treatment for enhanced user protection and has been developed from a sustainable perspective. That's why it has packaging made from FSC certified material and renewable sources. The OPGown comes with optimized packaging to ensure you get a great product in a smart dispenser using less material supporting your ambitions to reduce waste. evercare® SUPRA High Performance is reinforced with a breathable material in the critical zone, in the front and in the sleeves, to maximize protection against contamination of the surgical wound and to protect the user at wet procedures. All evercare® OP-Gowns provide an effective barrier to body fluids, skin particles and microbes and are designed for optimal comfort and ease of movement. The adjustable Velcro in the neck enables optimal fit while securing adjustments. High quality cuffs for improved comfort and tight fit round wrists. The gown is sterilized with EO-gas and packed in easy-to-open unit pouches. State-of-the-art packaging system with handy functional dispenser and labelling. Unit packs are supplied with barcoded removable labels for improved traceability.

Key Features

- ▶ Good protection with high breathability
- ▶ Super comfort with soft SMS
- ▶ Antistatic and alcohol repellent treated
- ▶ Very high abrasion resistance
- ▶ Extra soft stretchable neckline material for higher comfort
- ▶ Ultrasonic sealed sleeves to increase comfort by reducing excessive materials causing discomfort
- ▶ Gown with overlapping back and adjustable Velcro in the neck enables optimal fit for individual adjustments
- ▶ Clear stamp with size and performance type
- ▶ High quality cuffs for improved comfort and tight fit round wrists
- ▶ Belt card with clear indications of sterile / non-sterile for secure donning
- ▶ Book-folding enabling unpacking and donning without risk for contamination
- ▶ SMS material with low linting
- ▶ Comfortable design providing ease of movement
- ▶ All sizes available with or without 2 hand towels




REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
259011	SUPRA SP OP-Gown, sterile, Blue	Size M-short, Length 120 cm	21 / 42
259021	SUPRA SP OP-Gown, sterile, Blue	Size M, Length 131 cm	20 / 40
259031	SUPRA SP OP-Gown, sterile, Blue	Size L, Length 142 cm	20 / 40
259041	SUPRA SP OP-Gown, sterile, Blue	Size XL-short, Length 142 cm	20 / 40
259051	SUPRA SP OP-Gown, sterile, Blue	Size XL, Length 152 cm	18 / 36
259061	SUPRA SP OP-Gown, sterile, Blue	Size XXL, Length 158 cm	18 / 36
259071	SUPRA SP OP-Gown, sterile, Blue	Size XXL-Long, Length 170 cm	14 / 28
259081	SUPRA SP OP-Gown with hand towels, sterile, Blue	Size M-short, Length 120 cm	21 / 42
259091	SUPRA SP OP-Gown with hand towels, sterile, Blue	Size M, Length 131 cm	20 / 40
259101	SUPRA SP OP-Gown with hand towels, sterile, Blue	Size L, Length 142 cm	20 / 40
259111	SUPRA SP OP-Gown with hand towels, sterile, Blue	Size XL-short, Length 142 cm	20 / 40
259121	SUPRA SP OP-Gown with hand towels, sterile, Blue	Size XL, Length 152 cm	18 / 36
259131	SUPRA SP OP-Gown with hand towels, sterile, Blue	Size XXL, Length 158 cm	18 / 36
259141	SUPRA SP OP-Gown with hand towels, sterile, Blue	Size XXL-Long, Length 170 cm	14 / 28

Sales Unit: Transport Carton

Regulations

The product complies with legal requirements for medical devices. EO sterilized according to valid version of ISO 11135 standard.



REACH regulation:	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
Biological evaluation:	EN ISO 10993-7 EN ISO 10993-5 EN ISO 10993-10
Compliance with legal requirements:	Medical Device Regulation (EU) 2017/745, MDR Class Is
Compliance with product standards:	EN ISO 11135, EN 13795-1 Standard performance
Quality standards:	EN ISO 13485
Label information and packaging compliance:	EN ISO 20417, EN ISO 11607-1, EN ISO 15223-1

Material

Product Part	Raw Material & Characteristics
Body, sleeve, belt	PP SMS nonwoven, blue
Cuff	PES, rib knitted fabric
Neck ribbon	PP Spunbond nonwoven, white
Velcro band	PA woven ribbon with hook or loop
Belt card	Cellulose pulp (ECF), white
Wrapping crepe	Pulp (ECF), white
Sewing thread	PES, white
Hand towel	2pcs, Scrim Reinforced Tissue, white

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space below 25 °C. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 years

Packaging Information

Package	Material	Size (L x W x H) mm
Single Pack	Medical grade paper, PP/PA film	360 x 245 x 30
Dispenser Pack	Carton, TCF or ECF quality	381 x 279 x 320
Transport Carton	TCF or ECF quality, Carton	574 x 391 x 330
Packing Tape	PVC-free, Transparent	N/A

Barcodes / UDI-DI

REF	Single Pack	Dispenser Pack	Transport Carton
259011	06438129125344	06438129325348	06438129525342
259021	06438129125351	06438129325355	06438129525359
259031	06438129125368	06438129325362	06438129525366
259041	06438129125375	06438129325379	06438129525373
259051	06438129125382	06438129325386	06438129525380
259061	06438129125399	06438129325393	06438129525397
259071	06438129125405	06438129325409	06438129525403
259081	06438129125412	06438129325416	06438129525410
259091	06438129125429	06438129325423	06438129525427
259101	06438129125436	06438129325430	06438129525434
259111	06438129125443	06438129325447	06438129525441
259121	06438129125450	06438129325454	06438129525458
259131	06438129125467	06438129325461	06438129525465
259141	06438129125474	06438129325478	06438129525472

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.