



evercare SUPRA Special OP-Gowns, sterile

REF: 259151, 259161, 259171, 259181, 259191, 259201, 259211, 259221, 259231, 259241, 259251, 259261, 259271, 259281

evercare® SUPRA OP-Gown is a soft SMS gown designed to provide ultimate protection and super comfort. The gown is suitable for procedures with medium to high quantities of fluids and high comfort demands. SUPRA OP-Gown material has antistatic and alcohol repellent treatment for enhanced user protection and has been developed from a sustainable perspective. That's why it has packaging made from FSC certified material and renewable sources. The OP-Gown comes with optimized packaging to ensure you get a great product in a smart dispenser using less material supporting your ambitions to reduce waste. evercare® SUPRA High Performance is reinforced with a breathable material in the critical zone, in the front and in the sleeves, to maximize protection against contamination of the surgical wound and to protect the user at wet procedures. All evercare® OP-Gowns provide an effective barrier to body fluids, skin particles and microbes and are designed for optimal comfort and ease of movement. The adjustable Velcro in the neck enables optimal fit while securing adjustments. High quality cuffs for improved comfort and tight fit round wrists. The gown is sterilized with EO-gas and packed in easy-to-open unit pouches. State-of-the-art packaging system with handy functional dispenser and labelling. Unit packs are supplied with barcoded removable labels for improved traceability.

Key Features

- ▶ Good protection with high breathability
- ▶ Super comfort with soft SMS
- ▶ Antistatic and alcohol repellent treated
- ▶ Very high abrasion resistance
- ▶ Extra soft stretchable neckline material for higher comfort
- ▶ Reinforcement with breathable comfort layer on booth arms and chest
- ▶ Ultrasonic sealed sleeves to increase comfort by reducing excessive materials causing discomfort
- ▶ Gown with overlapping back and adjustable Velcro in the neck enables optimal fit for individual adjustments
- ▶ Clear stamp with size and performance type
- ▶ High quality cuffs for improved comfort and tight fit round wrists
- ▶ Belt card with clear indications of sterile / non-sterile for secure donning
- ▶ Book-folding enabling unpacking and donning without risk for contamination
- ▶ SMS material with low linting
- ▶ Comfortable design providing ease of movement
- ▶ All sizes available with or without 2 hand towels
- ▶ HP product variant is suitable for procedures with medium to high amount of fluids.




REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
259151	SUPRA HP OP-Gown, sterile, Blue	Size M-short, Length 120 cm	18 / 36
259161	SUPRA HP OP-Gown, sterile, Blue	Size M, Length 131 cm	16 / 32
259171	SUPRA HP OP-Gown, sterile, Blue	Size L, Length 142 cm	16 / 32
259181	SUPRA HP OP-Gown, sterile, Blue	Size XL-short, Length 142 cm	16 / 32
259191	SUPRA HP OP-Gown, sterile, Blue	Size XL, Length 152 cm	14 / 28
259201	SUPRA HP OP-Gown, sterile, Blue	Size XXL, Length 158 cm	14 / 28
259211	SUPRA HP OP-Gown, sterile, Blue	Size XXL-long, Length 170cm	10 / 20
259221	SUPRA HP OP-Gown with hand towels, sterile, Blue	Size M-short, Length 120 cm	18 / 36
259231	SUPRA HP OP-Gown with hand towels, sterile, Blue	Size M, Length 131 cm	16 / 32
259241	SUPRA HP OP-Gown with hand towels, sterile, Blue	Size L, Length 142 cm	16 / 32
259251	SUPRA HP OP-Gown with hand towels, sterile, Blue	Size XL-short, Length 142 cm	16 / 32
259261	SUPRA HP OP-Gown with hand towels, sterile, Blue	Size XL, Length 152 cm	14 / 28
259271	SUPRA HP OP-Gown with hand towels, sterile, Blue	Size XXL, Length 158 cm	14 / 28
259281	SUPRA HP OP-Gown with hand towels, sterile, Blue	Size XXL-long, Length 170 cm	10 / 20

Sales Unit: Transport Carton

Regulations

The product complies with legal requirements for medical devices. EO sterilized according to valid version of ISO 11135 standard.



REACH regulation:	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
Biological evaluation:	EN ISO 10993-7 EN ISO 10993-5 EN ISO 10993-10
Compliance with legal requirements:	Medical Device Regulation (EU) 2017/745, MDR Class Is
Compliance with product standards:	EN ISO 11135, EN 13795-1 High performance
Quality standards:	EN ISO 13485
Label information and packaging compliance:	EN ISO 20417, EN ISO 11607-1, EN ISO 15223-1

Material

Product Part	Raw Material & Characteristics
Body, sleeve, belt	PP SMS nonwoven, blue
Reinforcement sleeve and front	PP nonwoven / breathable PE-film, white
Cuff	PES, rib knitted fabric
Neck ribbon	PP Spunbond nonwoven, white
Velcro band	PA woven ribbon with hook or loop
Belt card	Cellulose pulp (ECF), white
Wrapping crepe	Pulp (ECF), white
Sewing thread	PES, white
Hand towel	2pcs, Scrim Reinforced Tissue, white

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space below 25 °C. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 years

Packaging Information

Package	Material	Size (L x W x H) mm
Single Pack	Medical grade paper, PP/PA film	360 x 245 x 30
Dispenser Pack	Carton, TCF or ECF quality	381 x 279 x 320
Transport Carton	TCF or ECF quality, Carton	574 x 391 x 330
Packing Tape	PVC-free, Transparent	N/A

Barcodes / UDI-DI

REF	Single Pack	Dispenser Pack	Transport Carton
259151	06438129125207	06438129325201	06438129525205
259161	06438129125214	06438129325218	06438129525212
259171	06438129125221	06438129325225	06438129525229
259181	06438129125238	06438129325232	06438129525236
259191	06438129125245	06438129325249	06438129525243
259201	06438129125252	06438129325256	06438129525250
259211	06438129125269	06438129325263	06438129525267
259221	06438129125276	06438129325270	06438129525274
259231	06438129125283	06438129325287	06438129525281
259241	06438129125290	06438129325294	06438129525298
259251	06438129125306	06438129325300	06438129525304
259261	06438129125313	06438129325317	06438129525311
259271	06438129125320	06438129325324	06438129525328
259281	06438129125337	06438129325331	06438129525335

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.