



evercare[®] Medical Face Mask type IIR
tieband,antifog,wide,white

REF: 2885

evercare[®] Medical face mask with tie-bands, type IIR provides a barrier to minimize the direct transmission of infective agents between medical staff and patient. It is a splash resistant mask offering protection against potentially contaminated fluid splashes. Recommended to be used during wet surgery.

Key Features

- ▶ EN 14683:2019
- ▶ Bacterial Filtration Efficiency (BFE) ≥ 98%
- ▶ Single use
- ▶ Flexible nose wire for proper fit
- ▶ Splash resistant
- ▶ Anti-fog
- ▶ Soft inner layer for comfort




REF	DESCRIPTION	SIZE	Pack Count dispenser/carton
2885	evercare [®] Medical Face Mask type IIR tie-band, antifog, wide, White	19 x 9 cm	50 / 300

Sales Unit: Transport Carton

Regulations

The product complies with legal requirements for medical devices.



REACH regulation:	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
Biological evaluation:	EN ISO 10993-5 EN ISO 10993-10
Compliance with legal requirements:	Medical Device Regulation (EU) 2017/745, MDR Class I
Compliance with product standards:	EN 14683:2019
Quality standards:	EN ISO 13485
Label information and packaging compliance:	EN ISO 15223-1, EN ISO 20417

Material

Product Part	Raw Material & Characteristics
1st Outer facing	Spunbond PP, white
2nd Outer facing	Spunbond PP, white
Filter media	Meltblown PP, white
Inner facing	Hydrophobic tissue, white
Tie-bands	Spunbond PP, white, left/right length 405 mm
Nose wire	PE coated metal wire, 135 mm
Antifog	Spunbond PP, white

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 years

Packing Information

Package	Material	Size (L x W x H) mm
Single Pack	N/A	N/A
Dispenser Pack	TCF or ECF quality	205 x 100 x 150
Transport Carton	TCF or ECF quality	320 x 225 x 320
Packing Tape	Transparent	N/A

Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
2885	N/A	06438129325072	06438129525076

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.