

evercare[®] Medical Face Mask type II ear-loop, green

REF: 2888

evercare[®] Medical face mask , standard with earloops, type II provides a barrier to minimize the direct transmission of infective agents between medical staff and patient. It is an allround face mask with high filtration capacity and low breathing resistance, designed to give an excellent /improved fit of the face mask for the user. Recommended to be used during dry surgery, with no risk of splashes.

Key Features

- ▶ EN 14683:2019
- ▶ Bacterial Filtration Efficiency (BFE) $\geq 98\%$
- ▶ Single use
- ▶ Flexible nose wire for proper fit
- ▶ Color: Outer facing green, inner facing white




Ear-loop

REF	DESCRIPTION	SIZE	Pack Count dispenser/carton
2888	evercare [®] Medical Face Mask type II ear-loop, Green	18 x 9 cm	50 / 300

Sales Unit: Transport Carton

Regulations

The product complies with legal requirements for medical devices.



REACH regulation:	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
Biological evaluation:	EN ISO 10993-5 EN ISO 10993-10
Compliance with legal requirements:	Medical Device Regulation (EU) 2017/745, MDR Class I
Compliance with product standards:	EN 14683:2019
Quality standards:	EN ISO 13485
Label information and packaging compliance:	EN ISO 15223-1, EN ISO 20417

Material

Product Part	Raw Material & Characteristics
Outer facing	Spunbond PP, green
Filter media	Meltblown PP, white
Inner facing	Thermalbonded PP, white
Ear loops	Nylon, Spandex & Polyester, 6mm wide, white
Nose wire	PE coated metal wire, 135 mm

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 years

Packing Information

Package	Material	Size (L x W x H) mm
Single Pack	N/A	N/A
Dispenser Pack	TCF or ECF quality	215 x 100 x 75
Transport Carton	TCF or ECF quality	235 x 220 x 245
Packing Tape	Transparent	N/A

Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
2888	N/A	06438129325195	06438129525199

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.