

KD-JECT® O

KD-JECT® O They are equipped with a connector incompatible to Luer connectors, thereby preventing accidental parenteral administration. They are intended for the enteral administration of liquid diets or drugs and for enteral aspiration or irrigation procedures. The scale allows exact measuring of the volumes to be administered or aspirated (accuracy of the scale in accordance with ISO 7886-1).


KD-JECT® O are non-toxic, pyrogen-free and have been sterilized by ethylene oxide (EO).



Product description

- Highly visible scale
- Siliconized for best gliding characteristics
- Large wings for easy and firm grip
- Piston stop that prevents the plunger being pushed out
- Purple plunger
- Color coded dispenser packaging which can be opened in two different ways
- PVC and DEHP-free, Latex free
- Shelf life 5 years

REF	Size (ml)	Volume scale (ml)	Nozzle	Tip	Pack count
876244	1	0,01	Centric	Oral tip	100
876251	3	0,1	Centric	Oral tip	100
876268	5	0,1	Centric	Oral tip	100
876275	10	0,2	Centric	Oral tip	100
876282	20	1	Centric	Oral tip	100



Component	Material	Packing material:	Sterilization:	Standards:
Barrel:	Ethylene Polypropylene Copolymer (PP)	Medical grade paper	Ethylene oxide (EO)	ISO 20695
Plunger:	Polyethylene (PE)	Transparent blister package PE/PA		ISO 7886-1
Rubber seal:	Synthetic rubber (polyisoprene)			CE 0123
Lubricant:	Silicone oil			class II a