

Embra® Optimia Nonwoven swabs, non-sterile

REF: 22300103, 22300203, 22300303, 22300403

Embra® Optimia Nonwoven swabs, 4-ply are produced from soft and highly absorptive nonwoven material. They are primarily used for cleaning and exudate absorption of wounds and abrasions.

Key Features

- ▶ 70% renewable materials (viscose)
- ▶ High quality viscose/polyester nonwoven material
- ▶ Low linting and high softness
- ▶ Keeps its shape, even when wet
- ▶ Free from binding agents and optical brighteners
- ▶ Effective in liquid and exudate absorption
- ▶ Less adherent to the wound than gauze devices
- ▶ Available in several sizes
- ▶ For external wound care, treating burns or abrasions, in the first aid area etc.
- ▶ Can be sterilised using EO
- ▶ Optimised packaging: convenient dispensing and better protection from contamination



REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
22300103	Optimia Nonwoven swabs, 4-ply, non-sterile, White	5x5 cm	150 / 13500
22300203	Optimia Nonwoven swabs, 4-ply, non-sterile, White	7,5x7,5 cm	150 / 9000
22300303	Optimia Nonwoven swabs, 4-ply, non-sterile, White	10x10 cm	150 / 5400
22300403	Optimia Nonwoven swabs, 4-ply, non-sterile, White	10x20 cm	75 / 900

Sales Unit: Transport Carton

Regulations

The product complies with legal requirements for medical devices.



REACH regulation:

The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).

Biological evaluation:

EN ISO 10993-5
EN ISO 10993-10

Compliance with legal requirements:

Medical Device Regulation (EU) 2017/745, MDR Class I

Compliance with product standards:

EN ISO 11737-1, EN 1644-1, EN 1644-2

Quality standards:

EN ISO 13485

Label information and packaging compliance:

EN ISO 15223-1, EN ISO 20417

Material

Product Part	Raw Material & Characteristics
Main material	Nonwoven, 70% viscose and 30% polyester, white

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space between 5 - 25 °C. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 years

Packaging Information

Package	Material	Size (L x W x H) mm
Single Pack	N/A	N/A
Dispenser Pack	TCF or ECF quality	22300103: 157 x 67 x 53 22300203: 232 x 67 x 77 22300303: 204 x 92 x 104 22300403: 204 x 92 x 104
Transport Carton	TCF or ECF quality	22300103: 490 x 280 x 420 22300203: 485 x 405 x 420 22300303: 430 x 330 x 570 22300403: 430 x 330 x 200
Packing Tape	PVC-free	N/A

Barcodes / UDI-DI

REF	Single Pack	Dispenser Pack	Transport Carton
22300103	N/A	06438129370218	06438129570212
22300203	N/A	06438129370225	06438129570229
22300303	N/A	06438129370232	06438129570236
22300403	N/A	06438129370249	06438129570243

Embra's mission is to develop responsible healthcare solutions focused on lowering the environmental impact while heightening quality for patients and caregivers. Together we can make a positive difference, so thank you for choosing Embra.