

flexGRIP®

catheter cushion and **stabilization device**

UNIVERSAL

the answer to known issues
with sutures:

- ✓ minimizes complications associated with suturing
- ✓ improves patient comfort and safety
- ✓ no more needle stick injuries by care givers



strong
comfortable
easy to use
versatile

flexGRIP®

catheter cushion and **stabilization device**

suitable for stabilizing & fixating a large range of catheters, such as:

nephrostomy (post-op),
drainage catheters, gripper
needles (Port-a-cath),
thrombolysis catheters,
jejunostomy, tubings, cables
and drivelines (eg. LVAD).



UNIVERSAL

features

- large diameter range: 5-18 Fr
- unmatched **strength** of securement
- **easy to use**; pre-mounted fixation leash
- additional lock for leash
- **skin-safe** hypoallergenic adhesive
- **soft** materials only

benefits

- universal system; wide variety of applications
- **prevent catheter migration** & related complications
- extra **safety** lip prevents accidental loosening
- **prevent MARSI** – medical adhesive related skin injury
- optimal **patient comfort**

ordering information

product No	description & size range	qty
UN511	FlexGRIP Universal 5-11 Fr	50 PCS/BOX
UN1118	FlexGRIP Universal 11-18 Fr	25 PCS/BOX

instruction video

