

# Embra<sup>®</sup> Surgical hood Gallant

REF: 277903

Embra<sup>®</sup> Headwears are designed to help prevent foreign particles such as hair or skin scales falling from medical personnel into the wound during surgical interventions. These very popular surgical hoods come with different special features and designs, such as integrated sweatbands for additional comfort, or open design for easier dressing placed on the head instead of over the head. Embra<sup>®</sup> surgical hood Gallant with shoulder protection is made in a soft and breathable viscose nonwoven material. Gallant integrates a sweatband for extra comfort and a tie-band for individual fit.

## Key Features

- ▶ Universal size
- ▶ Shoulder protection
- ▶ Integrated sweatband
- ▶ Tie-band
- ▶ 80% viscose, 20% polyester meaning most of the materials are from renewable resources
- ▶ Viscose from sustainable forests (FSC-certified)
- ▶ Production moved from Asia to Europe to further reduce emissions
- ▶ High-performance headwear protects from falling skin cells or stray hairs
- ▶ Breathable and comfortable headwear
- ▶ Drapes easily
- ▶ Suitable for long shifts in the most intense surgical environments
- ▶ Tested for skin irritation, skin sensitisation and in vitro cytotoxicity according to EN ISO 10993




REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
277903	Surgical hood Gallant, shoulder protection, tie-band, large sweatband, Green	Universal	50 / 300

Sales Unit: Transport Carton

## Regulations

The product complies with legal requirements for medical devices.



**REACH regulation:**

The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).

**Biological evaluation:**

EN ISO 10993-5  
EN ISO 10993-10

**Compliance with legal requirements:**

Medical Device Regulation (EU) 2017/745, MDR Class I

**Quality standards:**

EN ISO 13485

**Label information and packaging compliance:**

EN ISO 20417, EN ISO 15223-1

## Material

Product Part	Raw Material & Characteristics
Main material	Viscose nonwoven, green
Sweatband	Viscose spunlance nonwoven, white
Elastic band	Nylon and synthetic rubber, latex free, white
Tie bands	PP Spunbond, white
Sewing threads	Nylon and spun polyester, white

## Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 years

## Packaging Information

Package	Material	Size ( L x W x H ) mm
Single Pack	N/A	N/A
Dispenser Pack	TCF or ECF quality	185 x 300 x 145
Transport Carton	TCF or ECF quality	310 x 385 x 455
Packing Tape	Transparent, PVC-free	N/A

## Barcodes / UDI-DI

REF	Single Pack	Dispenser Pack	Transport Carton
277903	N/A	06438129370171	06438129570175

Embra's mission is to develop responsible healthcare solutions focused on lowering the environmental impact while heightening quality for patients and caregivers. Together we can make a positive difference, so thank you for choosing Embra.