

Article no. 72240
Category Protective sheet
Product Protective sheet light 80 x 175 cm
End use Effective protection of bed and bedlinen against body liquids



1. Ideal and hygienic as protection against body liquids.
2. Effective protection of bed and bedlinen.
3. Reduces the risk of contamination and improves hygiene.
4. Use as complement to other incontinence products, to prevent leaks.
5. Suitable for operation theatres, emergency care, ambulances etc. Also convenient to bring when travelling or on holiday.
6. Made of two layer tissue and one layer PE-film (polyethylene).
7. High-quality tissue, bleached. From 100% virgin fibers of TCF-pulp (Totally Chlorine Free).
8. Strong, resilient and liquid safe transparent PE-film.
9. Reinforced with 48 blue polyester threads.
10. With very high absorption abilities and strength.
11. The strongest sheet among our grades.
12. Using the right technique, the sheet is strong enough to slide a patient from one stretcher to another.
13. With leakage barriers along the sides, for higher security.
14. Size: 800 x 1750 mm.
15. CE-marked product. Fulfils the requirements of (EU) MDR 2017/745.
16. The sheets are packed in cartons on pallets. The pallets are stretched wrapped with transparent PE-film.
17. The cartons could be incinerated or recycled under standard conditions as cardboard.
18. The stretch film could be incinerated or recycled under standard conditions as standard plastic.
19. Produced in Sweden.

Size (cm)	Material	Units/ inner package Inner package/ transport unit	Transport units/ pallet Height of pallet (m)	Gross/Net weight carton (kg)	Measure of carton LxWxH (cm)	Pallet (EUR) size (mm)
80 x 175	Tissue + PE-film	100 1	24 1,090	7,60 7,23	48 29 23,5	1200 x 800

Finess Hygiene AB is certified in accordance with ISO 9001:2015, ISO 14001:2015 and ISO 13485:2016.

Our CE-marked products fulfils the requirements of (EU) MDR 2017/745.

2021-05-26 ver. 02