

# PRODUCT SHEET

## BODY WASH

1/2

### Soft pre-moistened bathing products

Swash<sup>®</sup> bathing products are especially designed to provide a daily rinse-free full body wash for care recipients with impaired mobility and the need of assisted bathing due to their health condition.

Swash<sup>®</sup> uses a system of wash gloves and wipes soaked in a special lotion which cleans and nourishes the skin. One package of 8 Swash<sup>®</sup> gloves and wipes is sufficient for one full body wash. Each section of the body is washed with a separate wash glove or wipe, then discarded with no need to rinse and/or dry the skin.



### Characteristics of the Swash<sup>®</sup> products:



• **Hypoallergenic lotion**, available with or without perfume



• **Free** from formaldehyde releasers



• **Dermatologically tested**



• Products can be **warmed up**



• **Nourishing ingredients** i.e. Vitamin E and Glycerin



• **pH skin-neutral** lotion



• **Nordic ecolabel** license number 5090 0291



• Made in **Holland** using a **fully automated production process**



• Production in a **clean room** environment guarantees high product quality



### Lyocell

The cloth of Swash<sup>®</sup> bathing products is, among other materials made from the wood fibre Lyocell. This natural fibre is characterised by its purity, softness and absorbency. Lyocell is produced in an environmentally friendly manner. The innovative production process that recycles both energy and water, ensures that the environment is burdened as little as possible.



## Product overview:

### Swash® bathing

Swash® products contain a balanced amount of lotion for optimal hygiene of the skin. The lotion thoroughly cleans and simultaneously nourishes the skin with caring ingredients, such as Vitamin E, Glycerin, Caprylyl Glycol and Glyceryl Caprylate. Swash® is produced with low-linting cloth. The combination of the fibres, Lyocell and polyester, ensures an optimal quality of the cloth.

### Wipes

Swash® Gold Wipes are made of a soft, medium weight cloth which provides a gentle and pleasant bath for people with vulnerable skin. The type of cloth used for the wipes is a cross-lapped Needle punch (non-woven).

### Gloves

Swash® gloves are made of firm, soft cloth based on cross-lapped Spunlace (non-woven) that has a light relief (embossing), which gives a pleasant feeling on the skin. The production process of the gloves uses an ultrasonic welding technique, which ensures a sufficiently strong glove. The glove makes it easier to clean difficult to reach body parts, such as skin folds and toes.

ITEM NUMBER	PRODUCT NAME	NUMBER OF PACKAGES PER BOX	NUMBER OF PACKAGES PER PALLET
A04028-8	<b>Swash® Gold Wipes</b> 8-pack Fragrance free	44	2.772
B04028-8	<b>Swash® Gold Wipes</b> 8-pack Fragranced	44	2.772
A04070-8	<b>Swash® Gold Gloves</b> 8-pack Fragrance free	40	2.520
B04070-8	<b>Swash® Gold Gloves</b> 8-pack Fragranced	40	2.520
A04070-6	<b>Swash® Gold Gloves</b> 6-pack Fragrance free	52	3.276

