

PRODUCT SHEET

INCONTINENCE CARE



1/2

Helps to prevent the occurrence of Incontinence-Associated Dermatitis (IAD)

Swash[®] Incontinence Care prevents the development of IAD and keeps it from spreading in its early stages, before the condition can worsen and pressure ulcers can occur. The extended line of Swash[®] Incontinence Care products can help the skin to recover faster thanks to its innovative and enhanced Trixtoric[®]-formula.

Swash[®] Perineum 3-in-1 incontinence care products are intended to **clean**, **care** for and **protect** the skin – in one simple step. **Swash[®] Perineum+** takes incontinence care one step further: The 4-in-1 incontinence care products are designed not only to **clean**, **care** for and **protect**, but also helps the skin to **restore**. This step is possible thanks to the newly developed lotion with Trixtoric[®] – a formula which helps the skin to restore faster.

The combination of several products into one simplifies incontinence care. While traditional incontinence care requires several different steps and products, thanks to Swash[®] Perineum+ ingredients, skin regimen can now be done in one single action. As a result, it is more comfortable for the care recipient, more convenient for the caregiver and encourages adherence to the care regimen.

Characteristics of Swash[®] incontinence products:



- Protective barrier with **3% Dimethicone**



- **100% fragrance-free** with neutral smell. Dermatologically **tested**.



- **pH skin-neutral** and **hypoallergenic** lotion



- Products can be **warmed**



- Contains **caring ingredients** such as Vitamin E, Hyaluronic Acid and Glycerine



- **Clean Room** production guarantees high quality product



- **Free** from formaldehyde releasers



- Made in **Holland** using a **fully automated production process**



- **Trixtoric[®]**: Enhanced formula exclusive to **Swash[®] Perineum+** helps the skin to restore



Lyocell

The cloth Swash[®] Perineum and Perineum+ products are, among other materials, made from the wood fibre Lyocell. This natural fibre is characterised by its purity, softness and absorbency. Lyocell is produced in an environmentally friendly manner. The innovative production process that recycles both energy and water, ensures that the environment is burdened as little as possible.

Product overview:

Swash® incontinence products

Swash® Perineum 3-in-1 products are moistened wash gloves or wipes infused with a skin **protectant**, 3% Dimethicone. This silicon-based oil leaves a thin layer on the skin, protecting it from humidity and irritation while being non-occlusive, hence allowing the skin to breath.

Swash® Perineum+ 4-in-1 products have the benefits of the 3-in-1 products plus the special **Trixotric®-Formula**: Swash® Perineum+ products put a protective, but breathable, layer on the skin, while supporting the skin's restoration with the help of the **Trixotric®-formula** – a combination of Hyaluronic Acid, Hydrogel and Natural Clay.

The perineal area can be cleaned thoroughly while being nourished with skin care ingredients such as Vitamin E, Glycerin, Hyaluronic acid, Capryl Glycol, Glyceryl Caprylate. Even severe contamination or incrustation can be removed with ease. Swash® is produced with a low-linting cloth and the combination of fibres, Lyocell and polyester, ensures the optimal quality of the products.

Swash® Cleansing Wipes are designed for a frequent cleansing of a soiled area. The lotion contains a soothing ingredient Jojoba which is gentle to the skin.

Swash® Perineum Wipes

The wipes are made of a soft, thick cloth based on a cross-lapped Needle punch (non-woven). The extra-long size of the cloth makes it possible to carry out incontinence clean-ups comfortably and quickly.

Swash® Perineum+ Wipes

The Wipes are made of a soft cloth Spunlace (non-woven) with a light relief for a good feel on the skin. Its special length helps to make incontinence hygienic, easy and comfortable.

Swash® Perineum Gloves

The gloves are made of a firm, soft cloth based on cross lapped Spunlace (non-woven) that has a light relief (embossing), which gives a pleasant feeling on the skin. The gloves are produced using an ultrasonic welding technique, which ensures a sufficiently strong glove is maintained. The glove makes it easier to clean difficult to reach body parts.

Swash® Cleansing Wipes

The wipes are made of a firm, soft cloth based on parallel Spunlace (non-woven).

ITEM NUMBER	PRODUCT NAME	NUMBER OF PACKAGES PER BOX / PALLET
Co4533-4	Swash® Perineum Wipes 4-pack (fragrance free)	64 / 4.032
Co4533-8	Swash® Perineum Wipes 8-pack (fragrance free)	40 / 2.520
Co4070-4	Swash® Perineum Gloves 4-pack (fragrance free)	80 / 5.040
Co4070-8	Swash® Perineum Gloves 8-pack (fragrance free)	40 / 2.520
Ko4540-4	Swash® Perineum+ Wipes 4-pack (fragrance free)	84 / 5.292
Ko4540-8	Swash® Perineum+ Wipes 8-pack (fragrance free)	56 / 3.528
Lo5012-48	Swash® Cleansing Wipes 48-pack (fragrance free)	12 / 756
Mo5012-48	Swash® Cleansing Wipes 48-pack (fragranced)	12 / 756
oo5010-50	Swash® Dry Bathing Gloves 50-pack	20 / 560

trixotric®