

PRODUCT SHEET

HAIR CARE

1/2



Rinse-free hair wash in bed!

Swash[®] hair care products are an inseparable component of Bathing 21[®]. The products are especially designed to provide a hygienic, comfortable hair wash for care recipients in need of assisted bathing. And that in 75% less time.

Compared to a traditional bed bath, Swash[®] hair care products make it easy to wash hair with no need to rinse. The products can be used every time the care recipient wants to wash his/her hair and as frequent as needed. It gives the care recipient a greater sense of self-esteem and self-confidence and it reduces the physical burden on the care recipient and the person providing care. These products further reduce the risk of contamination, waste and unnecessary spillage during washing.

Characteristics of Swash[®] hair care products:



With **fragrance** or **fragrance free** *



pH skin-neutral and **hypoallergenic** lotion



Dermatologically **tested**



Contains a **mild cleanser** i.e. Sodium Lauroyl Sarcosinate



Free from formaldehyde releasers



Products can be **warmed up**



Made in **Holland**



Production in a **clean room** environment guarantees high product quality



Lyocell

The cloth of Swash[®] Shampoo Cap is, among other materials made from the wood fiber Lyocell. This natural fiber is characterised by its purity, softness and absorbency. Lyocell is produced in an environmentally friendly manner. The innovative production process that recycles both energy and water, ensures that the environment is burdened as little as possible.



* Hypoallergenic perfume: none of the 26 fragrance ingredients that have been identified as allergens have been added.

Product overview:

Swash® Hair Care

The conventional hair washing, including preparation time, might last 20 till 25 minutes. With Swash® hair care products this takes only a few minutes because it reduces the need to turn or move the body. This ensures minimal stress for both care recipient and caregiver. The lotion thoroughly cleans the hair, is pH-skin neutral and has a positive effect on skin irritations.

Shampoo cap

Swash® Shampoo Cap is a single-use bathing cap, impregnated with a special hair care lotion. The cap has a diameter of 30 cm. The plasticised outside ensures that the lotion stays in the cap and protects the hands of the caregiver.

Hair lotion

Swash® Hair Lotion is a cleansing lotion, which contains nourishing ingredients. The bottle contains 100 ml lotion and has a disc-top cap that makes it easy to adapt the amount per hair wash. Especially suitable for care recipient with a head wound or to remove blood from the hair.



ITEM NUMBER	PRODUCT NAME	NUMBER OF PRODUCTS PER BOX	NUMBER OF PRODUCTS PER PALLET
Do4077-1	Swash® Shampoo Cap 1 cap Fragrance Free	52	3.276
Go4077-1	Swash® Shampoo Cap 1 cap Fragranced	52	3.276
Do4079-1	Swash® Hair Lotion 100 ml Fragrance Free	14	-