








HEMATURI KATETRAR

Hydrogelinkapslad latexkärna

Spets	Förstärkta		Nylon spiral	
	2-vägs	3-vägs	2-vägs	3-vägs
 Rund	–	–	30-50 ml Ch 18-24 2550H+Ch	30-50 ml Ch 18-24 2551H+Ch
 Couvelaire Whistle	30-50 ml Ch 18-24 1852H+Ch	30-50 ml Ch 18-24 1853H+Ch	30-50 ml Ch 18-24 2552H+Ch	30-50 ml Ch 18-24 2553H+Ch
 Dufour	30-50 ml Ch 18-24 1856H+Ch	30-50 ml Ch 18-24 1857H+Ch	30-50 ml Ch 18-24 2556H+Ch	30-50 ml Ch 18-24 2557H+Ch
 Dufour	60-80 ml Ch 18-24 1858H+Ch	60-80 ml Ch 18-24 1859H+Ch	–	–
 Delinotte	–	30-50 ml Ch 18-24 1861H+Ch	–	–
 * Le guillou	–	30-50 ml Ch 18-24 1864H+Ch	–	–
 * Urethral	–	30-50 ml Ch 18-24 1910VH+Ch	–	–

*Le Guillou och Urethral katetrarna har ett flertal dränageögon distalt om ballongen.