

PRODUCT DATA SHEET

TRANSOFIX (4090500)

Attribute Name	Attribute Value
Brand name	Transofix®
Product name (HC)	Transofix®
Product type	Transfer device
Application	Transferring fluids
Type of packaging	Blister
Packaging contains PVC	No
Packaging contains DEHP	No
Product (unpack.) contains PVC	No
Product (unpack.) cont. Latex	No
Product (unpack.) cont. DEHP	No
Packaging contains Latex	No
Product contains PVC	No
Product contains Latex	No
Product contains DEHP	No
Product contains Nickel	No

Attribute Name	Attribute Value
Product contains BPA	No
Sterility	sterile
Type of Sterilization	Gamma
Type of Sterilizat. (allowed)	not re-sterilizable
Disposable / Reusable	disposable
Priming volume	2.900 ml
Length	70.50 mm