

# PRODUCT DATA SHEET

INTROCAN SAFETY PUR 18G, 1.3X64-EU (4251620-01)



Attribute Name	Attribute Value
Brand name	Introcan Safety®
Product name (HC)	Introcan Safety®
Product type	IV catheter
Wing	No
Injection port	No
Flow rate	85.00 ml/min
Connection	Luer-Lock/female
Inner-ø catheter (metr.)	1.00 mm
Residual volume	0.140 ml
Color	green
Material of catheter	PUR
Safety clip	Yes

Attribute Name	Attribute Value
X-ray contrastable marking	X-ray contrast thread
Type of packaging	Blister
Product contains PVC	No
Product contains DEHP	No
Product contains Latex	No
Product contains Nickel	Yes
Remark critical attributes	Needle cont. Nickel
Disposable / Reusable	disposable
Catheter-Gauge	18
Outer-ø catheter (metr.)	1.30 mm
Length of catheter	64.0 mm
Length of catheter (inch)	2 1/2"