

PRODUCT DATA SHEET

Silicone mat (mesh), Standard 3/4, 402 x 250 mm, slip-free storage of instruments, good drying characteristics, gentle storage durin (JF940)



Attribute Name	Attribute Value
Product search type	Accessories Baskets
Aesculap Product Groups	Sterile Technology
Altern. Produktname	Positioning Mat
Type	silicone mat
Outer length (mm)	402 mm
Outer width (mm)	250 mm
Figure / Size	Standard 3/4
Surgical Step	Packaging
UDI Direct Marking	No
Disposable / Reusable	reusable
Product contains Latex	No
Sterility	non-sterile