

PRODUCT DATA SHEET

Bottom, without base perforation, Standard 3/4, outer length: 470 mm, outer width: 274 mm, outer height: 135 mm, inner length: 421 mm, inner width: 258 mm, inner height: 157 mm (JK742)



Attribute Name	Attribute Value
Product search type	Bottom
Aesculap Product Groups	Sterile Technology
Type	without base perforation
Lid Primeline	JP012, JP011, JP015, JP016, JP013, JP014
Lid Basis	JK785, JK786, JK789, JK787, JK788
Lid Vario	JK769
Outer length (mm)	470 mm
Outer width (mm)	274 mm
Outer height (mm)	135 mm
Inner length (mm)	421 mm
Inner width (mm)	258 mm
Inner height (mm)	157 mm
Outer lid optional Basis	JK775A, JK779A, JK776A, JK777A, JK778A
Outer Lid must-have Vario	JK755A, JK756A, JK757A, JK758A, JK759A
Figure / Size	Standard 3/4

Attribute Name	Attribute Value
Surgical Step	Packaging
Disposable / Reusable	reusable
Sterility	non-sterile
Product contains Latex	No