

PRODUCT SPECIFICATION

REF 1910 CEDERROTH 4-IN-1



1 GENERAL PRODUCT DESCRIPTION

Description

Cederroth 4-in-1 is a sterile universal dressing, which consists of two elastic bandages and one compress with a non-stick wound contact layer. The compress is stitched to one of the bandages. The dressing is wrapped in a pouch with printed pictorial user instructions.



The Cederroth 4-in-1 dressing is intended to be used in first aid situations. The dressing has four functions and thereby replaces a variety of traditional first aid dressings and compresses:

- Pressure pad for heavy bleeding
- Protective dressing for superficial wounds
- Burn dressing with special wound protecting surface
- Temporary support dressing

The product is intended for single use.

Storage information

Dry storage in room temperature.

For maximum protection, keep in original package until use.



2 DETAILED PRODUCT INFORMATION

Content per consumer unit

Product/Material

Big blood stopper in sterile package, packed in outer plastic pouch with illustrations

Quantity

1

Dimensions (w x l mm)

Bandage

3000 x 100

Compress

230 x 140

Material

Bandages: Viscose

Compresses: Absorbent layer: Viscose staple fibers

Non-adherent outer layer: Porous net made of high density polyethylene

Sterile barrier: Tyvek and transparent film made of polyethylene.

Outer pouch: Polyethylene with paper insert

Shelf life: 5 years from production date

3 LOGISTIC INFORMATION

	CU	SKU	Pallet
GTIN:	7310610019105	7310801910006	7312021910002
Content:	1 dressing	20 CU	14 SKU/layer x 6 layers
Size (w x l x h mm):	125 x 35x 65	205 x 260 x 175	800 x 1200 x 1200
Gross weight:	46 g	1100 g	117 kg
Packing material:	Plastic 6 g, Paper 5 g	Cardboard 156 g	-

