

# PRODUCT SPECIFICATION

REF 1911 CEDERROTH 4-IN-1 MINI



## 1 GENERAL PRODUCT DESCRIPTION

### Description

Cederroth 4-in-1 mini is a sterile universal dressing, which consists of one elastic bandage and one compress with a non-stick wound contact layer which is stitched to the bandage. The dressing is wrapped in a pouch with printed pictorial user instructions.

The Cederroth 4-in-1 mini dressing is intended to be used in first aid situations and it is especially suitable for fingers and toes. The dressing has four functions and thereby replaces a variety of traditional first aid dressings and compresses:

- Pressure pad for heavy bleeding
- Protective dressing for superficial wounds
- Burn dressing with special wound protecting surface
- Temporary support dressing

The product is intended for single use.

### Storage information

Dry storage in room temperature.

For maximum protection, keep in original package until use.



## 2 DETAILED PRODUCT INFORMATION

### Content per consumer unit

<u>Product/Material</u>	<u>Quantity</u>	<u>Dimensions (w x l mm)</u>	
		<u>Bandage</u>	<u>Compress</u>
Mini blood stopper in sterile package, packed in outer plastic pouch with illustrations	1	3000 x 60	120 x 80

### Material

Bandage:	Viscose
Compresses:	Absorbent layer: Viscose staple fibers Non-adherent outer layer: Porous net made of high density polyethylene
Sterile barrier:	Paper and transparent film made of polyethylene/polyamide.
Outer pouch:	Polyethylene with paper insert.

**Shelf life:** 5 years from production date

## 3 LOGISTIC INFORMATION

	CU	SKU	Pallet
<b>GTIN:</b>	7310610019112	7310801911003	7312021911009
<b>Content:</b>	1 dressing	30 CU	35 SKU/layer x 6 layers
<b>Size (w x l x h mm):</b>	40 x 25 x 65	155 x 210 x 105	800 x 1200 x 1080
<b>Gross weight:</b>	14 g	520 g	134 kg
<b>Packing material:</b>	Plastic 2 g, Paper 2 g	Cardboard 82 g	-



Version: 08  
Issued: 2018-02-27