

PRODUCT SPECIFICATION

REF 323700 SALVEQUICK WOUND CLEANSER 20-PACK



1 GENERAL PRODUCT DESCRIPTION

Description

Sterile sachets made of non-woven tissues, which are saturated in a wound cleansing liquid. The tissues are packed in a well-sealed foil, which prevents leakage and guarantees sterility. Salvequick Wound Cleanser Sachet is intended for mechanical cleansing of minor wounds when clean tap water is not available.

The product is a medical device, class IIa.

Storage information

Dry storage in room temperature. For maximum protection, keep in original package until use.



2 DETAILED PRODUCT INFORMATION

Content per consumer unit

<u>Product</u>	<u>Quantity</u>	<u>Dimensions (mm)</u>
Salvequick Wound Cleanser (Ref# 3239)	20	
Tissue		100 x 140
Sachet		60 x 58

Material

Tissue:	Nonwoven made of viscose
Liquid:	Physiological saline solution. Fill volume: 2,7 ml
Sachet:	Paper/LDPE/Aluminium/Surlyn®

Sterilization: Gamma irradiation

Shelf life: 5 years from production date

3 LOGISTIC INFORMATION

	CU:	SKU:	Pallet
GTIN:	7310613237001	7310803237002	-
Content:	20 sachets	20 CU	19 SKU / layer / 6 layers
Size (w x l x h mm):	62 x 66 x 77	135 x 330 x 160	Pallet 800 x 1200 x 960
Gross weight:	85 g	1841 g	236 kg
Packing material:	9 g	Carton 130 g	-



Orkla Care AB, Svetsarvägen 15, Box 1336, SE 171 26 Solna, Sweden, Phone +46 (0)10 142 64 00
www.firstaid.cederroth.com firstaid@cederroth.com

Version: 04
Issued: 2018-11-28



Teija Alander
SR QUALITY ASSURANCE & REGULATORY AFFAIRS MANAGER
WOUND CARE, ORKLA CARE AB

Teija alander
2018.11.28
12:31:23 +01'00'