

# PRODUCT SPECIFICATION

REF 51030127 SALVEQUICK DETECTABLE PLASTER



## 1 GENERAL PRODUCT DESCRIPTION

### Description

Detectable, adhesive bandages (plasters), of two sizes, medium and senior, packed in a cardboard refill intended for the Cederroth (Salvequick) dispenser for adhesive bandages.

The Salvequick Detectable bandages are detectable by most common metal detectors and the blue color also makes them easy to detect visually, which makes these plasters suitable for restaurants and the food industry. The plasters can easily be pulled out with just one hand and applied on the wound.

The adhesive bandages consist of a backing, evenly spread with an adhesive coating and with an absorbent pad fixed centrally to the dressing. When the adhesive bandage is pulled out from the package, the lower cover will remain behind. Thus one half of the bandage can be applied with one hand and after removing the rest of the cover the plaster is ready for use.

### Storage information

Dry storage in room temperature. For maximum protection, keep in original package until use.



## 2 DETAILED PRODUCT INFORMATION

### Content per consumer unit

<u>Product/Material</u>	<u>Quantity (pcs)</u>	<u>Dimensions (mm)</u>	
Refill packages	6	30 x 88 x 100	
<u>Per refill package:</u>		<u>Dressing</u>	<u>Wound Pad</u>
Medium plasters	21	72 x 19	24 x 12
Senior plasters	14	72 x 25	24 x 17

### Material

Film/backing:	Blue, soft and flexible polyurethane/polystyrene foil
Adhesive:	Acrylate
Wound pad:	Laminate of viscose, polypropylene, polyethylene with polyethylene wound contact layer and a metal detectable aluminium foil
Facing:	Siliconized paper and plastic
Protective paper:	Cold seal paper coated with latex-based adhesive
Wrapping:	Carton
<b>Shelf life:</b>	5 years

## 3 LOGISTIC INFORMATION

	CU	SKU	Pallet
<b>GTIN:</b>	7310610019327	7310800019403	7312020020085
<b>Content:</b>	6 refills	10 CU	11 SKU/layer x 5 layers
<b>Size (w x l x h mm):</b>	125 x 100 x 90	130 x 513 x 190	Pallet 800 x 1200 x 1100
<b>Gross weight:</b>	188 g	2093 g	140 kg
<b>Packing material:</b>	Carton 25 g	Corrugated board 213 g	-

