PRODUCT SPECIFICATION

REF6454SALVEQUICK TEXTILE (FABRIC) FINGERTIP XL
PLASTERS (ADHESIVE BANDAGES)

1 GENERAL PRODUCT DESCRIPTION

Description

Adhesive, butterfly-shaped, fingertip bandages (plasters) packed in a cardboard refill package, intended for the Cederroth (Salvequick) dispenser for adhesive bandages.

The plasters can easily be pulled out with just one hand and applied on the wound.

The adhesive bandage consists of a backing, evenly spread with an adhesive coating. An absorbent pad is fixed centrally to the dressing. The adhesive is covered with a release-coated foil that consists of two parts partly overlapping each other. When the adhesive bandage is taken out from the package, the release-coated plastic foil remains behind. One half of the dressing can now be applied on the skin using one hand only. The rest of the coated paper can then be removed and the adhesive bandage is ready for use.

Storage information

Dry storage in room temperature. For maximum protection, keep in original package until use.

2 DETAILED PRODUCT INFORMATION

Content per consumer unit

Product/Material	<u>Quantity (pcs)</u>	<u>Dimensions (</u>	<u>(mm)</u>
Refill packages	6	30 x 88 x100	
<u>Per refill package:</u>	15	<u>Dressing</u>	<u>Wound Pad</u>
Fingertip bandage		87 x 56 (23)	16 x 39
Material Film/backing: Adhesive Wound pad Facing Protective paper Wrapping	Elastic, beige viscose/polyamide fabric with hydrophobic treatment Acrylate Viscose, polypropylene, polyethylene with polyethylene wound contact layer Siliconized paper and plastic Coldseal paper coated with latex-based adhesive Carton		

3 LOGISTIC INFORMATION

Shelf life:

	CU	SKU	Pallet
GTIN:	7310610064549	7310806454000	7312020064546
Content:	6 refills	10 CU	11 SKU/layer x 5 layers
Size (w x l x h mm):	150 x 100 x 90	155 x 513 x 190	Pallet 800 x 1200 x 1100
Gross weight:	194 g	2177 g	145g
Packing material:	Carton 28 g	Corrugated board 237 g	-



Version: 05 Issued: 2015-09-29







5 years