

# Curan Lady

Compact ready to use female device



The Curan Lady is unique by its compact, attractive and elegant design. The catheters are ready to use and extremely user friendly. The catheters are supplied in individual, sterile, single-use packages and free from the aggressive softening agent DEHP. Each catheter is individually tested and guaranteed leakage-free. All of our catheters feature polished eyelets, making insertion risk-free and preventing irritation. All urine collection bags available on the market can be used in combination with this product.



- 30 pieces packaged



- Polished Eyelets



- Compact Design

## The main advantages at a glance

- Free from the aggressive softening agent DEHP
- Smooth, polished eyelets for irritation-free insertion
- Guaranteed leakage-free
- Attractive, compact and elegant design
- Ready to use
- Extremely user friendly

## Female catheters

The female catheters are 3.35 inches/8.5 cm long and are available in the following charriere sizes:

### Product code

CL 08 Compact ready to use female device  
CL 10 Compact ready to use female device  
CL 12 Compact ready to use female device  
CL 14 Compact ready to use female device

### Size

Ch 08  
Ch 10  
Ch 12  
Ch 14

### Pack

box of 30  
box of 30  
box of 30  
box of 30

### Colour code

● Blue  
● Black  
● White  
● Green

Distributor:

**Curan: specialist in intermittent catheterization and urology products**

