

Curan Lady

Compact ready to use female device

lady



The Curan Lady is unique by its compact, attractive and elegant design. The catheters are ready to use and extremely user friendly. The catheters are supplied in individual, sterile, single-use packages and free from the aggressive softening agent DEHP. Each catheter is individually tested and guaranteed leakage-free. All of our catheters feature polished eyelets, making insertion risk-free and preventing irritation. All urine collection bags available on the market can be used in combination with this product.



- 30 pieces packaged



- Polished Eyelets



- Compact Design

The main advantages at a glance

- Free from the aggressive softening agent DEHP
- Smooth, polished eyelets for irritation-free insertion
- Guaranteed leakage-free
- Attractive, compact and elegant design
- Ready to use
- Extremely user friendly

Female catheters

The female catheters are 3.35 inches/8.5 cm long and are available in the following carrier sizes:

Product code

CL 08 Compact ready to use female device
CL 10 Compact ready to use female device
CL 12 Compact ready to use female device
CL 14 Compact ready to use female device





Size

Ch 08
Ch 10
Ch 12
Ch 14

Pack

box of 30
box of 30
box of 30
box of 30

Colour code

 Blue
 Black
 White
 Green

Distributor:

Curan: specialist in intermittent catheterization and urology products

