

# Curan Catheters Plain

For intermittent self-catheterization. DEHP-free



The Curan catheters plain are smooth, uncoated and free from the aggressive softening agent DEHP. They are packaged individually and sterile, for single use. Both the male and the female catheters are fitted out with smooth polished eyelets, which makes insertion risk-free and prevents irritation.



• Colour coding



• Polished Eyelets



• Flexible connector

## The main advantages at a glance

- Free from the aggressive softening agent DEHP
- Made according to the wishes of users and professionals
- Based on years of knowledge and experience
- Extensively tested, so safe and reliable

## Female catheters

The female catheters are 6 inches long and are available in the following charriere sizes:

Product code	Size	Pack	Colour code
F 8 Female plain Straight tip	Ch 8	box of 30	● Blue
F 10 Female plain Straight tip	Ch 10	box of 30	● Black
F 12 Female plain Straight tip	Ch 12	box of 30	● White
F 14 Female plain Straight tip	Ch 14	box of 30	● Green
F 16 Female plain Straight tip	Ch 16	box of 30	● Orange
F 18 Female plain Straight tip	Ch 18	box of 30	● Red

## Male catheters

The male catheters are 16 inches long and are available in the following charriere sizes:

Product code	Size	Pack	Colour code
M 8 Male plain Straight tip	Ch 8	box of 30	● Blue
M 10 Male plain Straight tip	Ch 10	box of 30	● Black
M 12 Male plain Straight tip	Ch 12	box of 30	● White
M 12C Male plain Tiemann tip	Ch 12	box of 30	● White
M 14 Male plain Straight tip	Ch 14	box of 30	● Green
M 14C Male plain Tiemann tip	Ch 14	box of 30	● Green
M 16 Male plain Straight tip	Ch 16	box of 30	● Orange
M 18 Male plain Straight tip	Ch 18	box of 30	● Red

Distributor:

**Curan: specialist in intermittent catheterization and urology products**

