

Curan Advantage Plus Hydrophilic Coated Catheters

Ready to use with water pocket



The Curan Advantage Plus hydrophilic coated catheters are ready to use thanks to the extra element: the water pocket. This leak-free water pocket, filled with sterile water, can be activated with a light squeeze. The sterile water activates the hydrophilic coating within 25 seconds. This promotes additional comfort during insertion. Thanks to the handy hanger sticker, the packaging can be placed vertically in any room during activation of the comfort coating. In the meantime, the user can prepare for catheterization. Furthermore, the male catheters are fitted with our unique Blue Grip®; an ingenious tool for inserting the catheter in a hygienic way, without having to touch it with the hands. The Blue Grip® is easy, more hygienic and significantly decreases the risk of infection. The Curan Advantage Plus is free from the aggressive softening agent DEHP. They are packaged individually and sterile, for single use. The Curan Catheters are fitted with smoothly polished eyelets, making insertion risk-free and preventing irritation.



• Water pocket



• Blue Grip®



• Polished eyelets

The main advantages at a glance

- Free from the aggressive softening agent DEHP
- Ready to use thanks to the additional water pocket
- The water pocket can be activated with a simple squeeze
- The sterile water activates the comfort coating
- Guaranteed leakage-free
- Handy hanger sticker for vertical placement in any room during activation time
- Easy and hygienic insertion (less risk of infection) thanks to the unique Blue Grip®
- Fitted out with smoothly polished eyelets, which prevents irritation during insertion

Female catheters

The female catheter are 6 inches long and are available in the following charriere sizes:

| Product code | Size | Pack | Colour code |
|--------------------------------|-------|-----------|-------------|
| WHF 8 Curan coated RTU Female | Ch 8 | box of 25 | ● Blue |
| WHF 10 Curan coated RTU Female | Ch 10 | box of 25 | ● Black |
| WHF 12 Curan coated RTU Female | Ch 12 | box of 25 | ● White |
| WHF 14 Curan coated RTU Female | Ch 14 | box of 25 | ● Green |
| WHF 16 Curan coated RTU Female | Ch 16 | box of 25 | ● Orange |

Male catheters

The male catheters, including our unique Blue Grip®, are 16 inches long and are available in the following charriere sizes:

| Product code | Size | Pack | Colour code |
|--|-------|-----------|-------------|
| WHM 8 Curan coated RTU + Bluegrip® Male | Ch 8 | box of 25 | ● Blue |
| WHM 10 Curan coated RTU + Bluegrip® Male | Ch 10 | box of 25 | ● Black |
| WHM 10C Curan coated RTU+Bluegrip® & Tiemann tip | Ch 10 | box of 25 | ● Black |
| WHM 12 Curan coated RTU + Bluegrip® Male | Ch 12 | box of 25 | ● White |
| WHM 12C Curan coated RTU+Bluegrip® & Tiemann tip | Ch 12 | box of 25 | ● White |
| WHM 14 Curan coated RTU + Bluegrip® Male | Ch 14 | box of 25 | ● Green |
| WHM 14C Curan coated RTU+Bluegrip® & Tiemann tip | Ch 14 | box of 25 | ● Green |
| WHM 16 Curan coated RTU + Bluegrip® Male | Ch 16 | box of 25 | ● Orange |

Distributor:

Curan: specialist in intermittent catheterization and urology products

