



# 87111

## TOE kit, Polyurethane (incl. biteguard, syringe and transmission gel)

SIZE: **3 x 100 cm**

QUANTITY: **12**

CONTAINS LATEX: **No**

STERILE: **Yes**

STERILIZATION METHOD: **GAMMA**

MATERIAL: **Polyurethane**

GTIN PRIMARY: **50748426120463**

GTIN SALES: **50748426120463**

COMMODITY CODE: **4014100000**

INTENDED USE: **to cover an ultrasound transducer, used in natural body orifices, less than 1 hour**

SHELF LIFE: **3 Years**

DEHP/PHTALATES: **No**

GMDN CODE: **60608**

CE NOTIFYING BODY: **CE0482**

MEDICAL CLASSIFICATION: **Ila**

BOX LENGTH (MM): **380**

BOX WIDTH (MM): **290**

BOX HEIGHT: **140**

BOX M3: **0.015428**

BOX WEIGHT (GRAM): **950**

LABEL LANGUAGES: **DA,DE,ES,FI,FR,IT,NL,PL,PT,SE**

---

2X STERILISATION ALLOWED?: **No**

---

ADDRESS OF LEGAL MANUFACTURER: **Hekkehorst 24, 7207 BN Zutphen, The Netherlands**

---

BISPHENOL A: **No**

---

IFU: **Yes**

---

SINGLE USE: **Yes**

---

SPECIAL STORAGE CONDITIONS?: **Keep Away From Sunlight Symbol**

---

CASES PER PALLET CONVERSION: **120**

---

INFORMATION LAST UPDATED ON: **21 Jun 2022**

---