

**88662****PIP Ultracover™**SIZE: **20 x 300 mm**QUANTITY: **100**CONTAINS LATEX: **No**STERILE: **No**STERILIZATION METHOD: **Non Sterile**MATERIAL: **Polyisoprene**FOLDED/ROLLED: **Rolled**GTIN PRIMARY: **50748426100465**GTIN SALES: **50748426100465**COMMODITY CODE: **4014100000**INTENDED USE: **to cover an ultrasound transducer, used in natural body orifices, less than 1 hour**SHELF LIFE: **3 Years**DEHP/PHTHALATES: **No**GMDN CODE: **60609**CE NOTIFYING BODY: **CE0482**MEDICAL CLASSIFICATION: **Ila**BOX LENGTH (MM): **290**BOX WIDTH (MM): **195**BOX HEIGHT: **120**BOX M3: **0.006786**

BOX WEIGHT (GRAM): **800**

---

LABEL LANGUAGES: **BG,CZ,DE,DK,EE,EN,ES,FI,FR,GR,HR,HU,IS,IT,LT,LV,NL,NO,PL,PT,RO,RU,SE,SI,SK,SR,TR**

---

2X STERILISATION ALLOWED?: **No**

---

ADDRESS OF LEGAL MANUFACTURER: **Hekkehorst 24, 7207 BN Zutphen, The Netherlands**

---

BISPHENOL A: **No**

---

IFU: **Yes**

---

SINGLE USE: **Yes**

---

SPECIAL STORAGE CONDITIONS?: **Keep Away From Sunlight Symbol**

---

CASES PER PALLET CONVERSION: **240**

---

INFORMATION LAST UPDATED ON: **21 Jun 2022**

---