




Intersurgical InTube™ tracheal tubes


A comprehensive range of ET tubes for adult and paediatric use



Airway Management • Airway Devices

Intersurgical InTube™ tracheal tubes

Designed for use in anaesthesia, A&E and emergency medicine, the Intersurgical InTube range offers three varieties of tracheal tube: cuffed, reinforced and uncuffed, in a comprehensive range of sizes from paediatric to adult. Supplied sterile, individually packed, latex free and for single use.





illustrated

Cuffed

Code	Description	Size (Inner Diameter)	Box Qty.	
8040100	InTube tracheal tube, cuffed	10.0mm	10	
8040095	InTube tracheal tube, cuffed	9.5mm	10	
8040090	InTube tracheal tube, cuffed	9.0mm	10	
8040085	InTube tracheal tube, cuffed	8.5mm	10	
8040080	InTube tracheal tube, cuffed	8.0mm	10	
8040075	InTube tracheal tube, cuffed	7.5mm	10	
8040070	InTube tracheal tube, cuffed	7.0mm	10	
8040065	InTube tracheal tube, cuffed	6.5mm	10	
8040060	InTube tracheal tube, cuffed	6.0mm	10	
8040055	InTube tracheal tube, cuffed	5.5mm	10	
8040050	InTube tracheal tube, cuffed	5.0mm	10	
8040045	InTube tracheal tube, cuffed	4.5mm	10	

Make an enquiry



illustrated

Wire-reinforced

Code	Description	Size (Inner Diameter)	Box Qty.	
8060100	InTube tracheal tube, cuffed wire-reinforced	10.0mm	10	
8060095	InTube tracheal tube, cuffed wire-reinforced	9.5mm	10	
8060090	InTube tracheal tube, cuffed wire-reinforced	9.0mm	10	
8060085	InTube tracheal tube, cuffed wire-reinforced	8.5mm	10	
8060080	InTube tracheal tube, cuffed wire-reinforced	8.0mm	10	
8060075	InTube tracheal tube, cuffed wire-reinforced	7.5mm	10	
8060070	InTube tracheal tube, cuffed wire-reinforced	7.0mm	10	
8060065	InTube tracheal tube, cuffed wire-reinforced	6.5mm	10	
8060060	InTube tracheal tube, cuffed wire-reinforced	6.0mm	10	
8060055	InTube tracheal tube, cuffed wire-reinforced	5.5mm	10	
8060050	InTube tracheal tube, cuffed wire-reinforced	5.0mm	10	
8060045	InTube tracheal tube, cuffed wire-reinforced	4.5mm	10	

Make an enquiry



illustrated

Uncuffed

Code	Description	Size (Inner Diameter)	Box Qty.
8050070	InTube tracheal tube, uncuffed	7.0mm	10
8050065	InTube tracheal tube, uncuffed	6.5mm	10
8050060	InTube tracheal tube, uncuffed	6.0mm	10
8050055	InTube tracheal tube, uncuffed	5.5mm	10
8050050	InTube tracheal tube, uncuffed	5.0mm	10
8050045	InTube tracheal tube, uncuffed	4.5mm	10
8050040	InTube tracheal tube, uncuffed	4.0mm	10
8050035	InTube tracheal tube, uncuffed	3.5mm	10
8050030	InTube tracheal tube, uncuffed	3.0mm	10
8050025	InTube tracheal tube, uncuffed	2.5mm	10

Make an enquiry

Sterile

For the full range of Airway Management devices and accessories visit

www.intersurgical.com/products/airway-management



Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK
 T: +44 (0)118 965 6300 F: +44 (0)118 965 6356 info@intersurgical.com www.intersurgical.com



The manufacturer Intersurgical Ltd is certified to ISO 9001:2015, ISO 13485:2016 and ISO 14001:2015

Please think before you print
 Save energy and paper.
 If you must print this information sheet please print it double sided.

IS6.4 • Issue 7 07.19

**UK ▪ Ireland ▪ France ▪ Germany ▪ Spain ▪ Portugal ▪ Italy ▪ Benelux ▪ Sweden ▪ Denmark ▪ Lithuania ▪ Russia ▪ Czech Republic
 Turkey ▪ South Africa ▪ China ▪ Japan ▪ Taiwan ▪ Philippines ▪ USA ▪ Canada ▪ Colombia ▪ Australia**