

Product Name	i-gel adult, supraglottic airway		
Reference numbers	8203000	8204000	8205000
Product Type	i-gel Supraglottic Airway		
Manufacturer	Intersurgical		
Materials	PP - Protective cradle TPE - Body ABS - Connector		
Packaging			
Unit Packaging	Individually packaged		
Box quantity	25		
Pallet quantity	800		
Instructions for Use	Instructions for use with product Instructions for use with each box NO YES		
Certification	CE0120	Class 2A	
Applicable Standards	BS EN ISO 13485: Medical devices — Quality management systems — Requirements for regulatory purposes ISO 9001: Quality Management Systems - Requirements ISO 15223-1: Medical Devices - Symbols to be Used with Medical Device Labels, Labelling, & Information to be Supplied ISO 11712: Anaesthetic and Respiratory equipment - Supralaryngeal airways and connectors		
Sterile	Yes		
Shelf life	3 Years		
Disposal Methods	As hospital guidelines		
Re-use	 Single use only		

Product characteristics

Technical information	8203000	8204000	8205000
Patient	Small	Medium	Large
Size	3	4	5
Patient weight (kg)	30 - 60	50 - 90	90+
Weight - including packaging (g)	102	127	149
Internal volume of airway channel (ml)	23	27	34



**Material
Abbreviations**

ABS	Acrylonitrile Butadiene Styrene
AL	Aluminium
HDPE	High Density Polyethylene
HIPS	High Impact Polystyrene
LDPE	Low Density Polyethylene
MBS	Methacrylate Butadiene Styrene
NRL	Natural Rubber Latex
PC	Polycarbonate
PE	Polyethylene
PP	Polypropylene
PVC	Poly Vinyl Chloride
SB	Styrene Butadiene
Si	Silicone
SS	Stainless Steel
TPE	Thermoplastic Elastomer
SBC	Styrene Butadiene Copolymer

