

Product Name	i-gel paediatric, supraglottic airway		
Reference numbers	8201000	8215000	8202000 8225000
Product Type	i-gel Supraglottic Airway		
Manufacturer	Intersurgical		
Materials	PP - Protective cradle TPE - Body ABS - Connector		
Packaging	Unit Packaging Individually packaged		
Box quantity	10		
Pallet quantity	1020 (8225000 - 800)		
Instructions for Use	Instructions for use with product	NO	
	Instructions for use with each box	YES	
Certification	CE0120	Class 2A	
Applicable Standards	BS EN ISO 13485: Medical devices — Quality management systems — Requirements for regulatory purposes		
	ISO 9001: Quality Management Systems - Requirements		
	ISO 15223-1: Medical Devices - Symbols to be Used with Medical Device Labels, Labelling, & Information to be Supplied		
	ISO 11712: Anaesthetic and Respiratory equipment - Supralaryngeal airways and connectors		
Sterile	Yes		
Shelf life	2 years		
Disposal Methods	As hospital guidelines		
Re-use	Single use only		

Product characteristics

Technical information

	8201000	8215000	8202000	8225000
Patient size	Neonate	Infant	Small Paediatric	Large Paediatric
i-gel size	1	1.5	2	2.5
Patient weight (kg)	2-5	5-12	10-25	25-35
Weight - including packaging (g)	38	52	65	91
Internal volume of airway channel (ml)	6	8	14	21



**Material
Abbreviations**

ABS	Acrylonitrile Butadiene Styrene
AL	Aluminium
HDPE	High Density Polyethylene
LDPE	Low Density Polyethylene
PC	Polycarbonate
PE	Polyethylene
PP	Polypropylene
PVC	Poly Vinyl Chloride
SB	Styrene Butadiene
Si	Silicone
SS	Stainless Steel
TPE	Thermoplastic Elastomer

