

Product Name	i-gel O2 Resus Pack		
Reference numbers	8703000	8704000	8705000
Product Type	i-gel® supraglottic airway		
Manufacturer	Intersurgical		
Materials	TPE - I-gel O2 body Lubricating jelly - OptiLube ABS - I-gel O2 connector PU foam - Airway support strap PP - Protective cradle PVC - Naso-gastric tube		
Packaging			
Unit Packaging	Individually packaged		
Box quantity	6		
Pallet quantity	612		
Instructions for Use	Instructions for use with product NO Instructions for use with each box YES		
Certification	CE0120	Class 2A	
Applicable Standards	BS EN ISO 13485: Medical devices — Quality management systems — Requirements for regulatory purposes ISO 9001: Quality Management Systems - Requirements ISO 15223-1: Medical Devices - Symbols to be Used with Medical Device Labels, Labelling, & Information to be Supplied ISO 11712: Anaesthetic and Respiratory equipment - Supralaryngeal airways and connectors BS EN ISO 8836: Suction Catheters for Use in the Respiratory Tract		
Sterile	Yes, supplied sterile		
Shelf life	2 years		
Disposal Methods	As hospital guidelines		
Re-use	② Single use only		


Product characteristics

Technical information

	8703000	8704000	8705000
i-gel O2 size	3	4	5
Recommended patient weight (kg)	30-60	50-90	90 +
Hook ring colour code	Yellow	Green	Orange
Complete I-gel O2 resus pack weight (g)	129	146	170

Technical information - Naso-gastric tube

Shaft length (mm)	1200 +/- 40
Outside diameter (mm)	3.95 +/- 0.10
Inside diameter (mm)	2.6 +/- 0.15



*Width and Length correspond to i-gel O2 not packet size

The i-gel O2 resus pack contains:

- i-gel O2 supraglottic airway
- OptiLube lubricant sachet
- Naso-gastric tube
- Airway support strap



**Material
Abbreviations**

ABS	Acrylonitrile Butadiene Styrene
AL	Aluminium
HDPE	High Density Polyethylene
LDPE	Low Density Polyethylene
PC	Polycarbonate
PE	Polyethylene
PP	Polypropylene
PVC	Poly Vinyl Chloride
SB	Styrene Butadiene
Si	Silicone
SS	Stainless Steel
TPE	Thermoplastic Elastomer

