



## Technical data sheet

**Product category:** Sterile Syringes with NaCl 0,9%  
**Description:** Syringe with NaCl 0,9%, luer lock for flushing only



**Technical drawing:**



**Material:** Syringe made of polypropylene  
 Syringe cap made of styrene block copolymer  
 Sodium chloride 0,9%

**Sterilization:** Radiation

**Brand:** Sterisets Medical Products

**Product overview:** Sterile single pack

REF	Description	Quantity in pcs	
		Inner box	Transport carton
14361S	Syringe Luer Lock 3 ml NaCl 0,9%, for flushing only *S	60	480
14362S	Syringe Luer Lock 5 ml NaCl 0,9%, for flushing only *S	60	480
14363S	Syringe Luer Lock 10 ml NaCl 0,9%, for flushing only *S	40	240
14364S	Syringe Luer Lock 20 ml NaCl 0,9%, for flushing only *S	25	150