

SELEFA® Apron with sleeves

REF: 100028

SELEFA® Apron with sleeves is an impermeable protective disposable gown with thumb loops made of liquid-resistant polyethylene. It is a PPE Cat I product, designed to protect clothing from spills and fluid splashes. Light weight and flexible design allows for comfortable wear.

Key Features

- ▶ Open back
- ▶ Thumb loops at cuffs to prevent sleeves from riding up when donning gloves
- ▶ Quick to put on and easy to take off by splitting at the neck
- ▶ Individually folded and individually bagged




REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
100028	SELEFA® Apron with sleeves, Light blue	110 x 135 cm	20 / 160

Sales Unit: Transport Carton

Regulations

The product complies with legal requirements for personal protective equipments.



REACH regulation:	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
Compliance with legal requirements:	Personal Protective Equipment Regulation (EU) 2016/425, PPE R category I. EU Declaration of Conformity for Personal Protective Equipment: https://www.evercare-medical.com/about-us
Label information and packaging compliance:	EN ISO 15223-1

Material

Product Part	Raw Material & Characteristics
Main material	LLDPE with LDPE, 0,025 mm

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 years

Packaging Information

Package	Material	Size (L x W x H) mm
Single Pack	PE	
Dispenser Pack	Corrugated cardboard	200 x 170 x 160
Transport Carton	Corrugated cardboard	695 x 410 x 175
Packing Tape	PVC-free, Transparent	N/A

Barcodes / UDI-DI

REF	Single Pack	Dispenser Pack	Transport Carton
100028	N/A	6438129800135	6438129800142

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.