

## evercare<sup>®</sup> inLine 3-way stopcock, Extension line, PVC-free

REF: 10548-01, 10549-01, 10550-01, 10551-01, 10552-01, 10553-01, 10713

evercare<sup>®</sup> inLine 3-way stopcocks are a and broad range of product for different needs. 3-way stopcock with extension tubing in different legths enabling connections to two sources of infusion and/or injection or a combination of both. Stopcock are sterile, single use products that are used to conduct intravenous fluids as an extension of the infusion or perfusion device, by enabling connections to two sources of infusion and an infusion device with vascular access.

### Key Features

- ▶ Lipid resistant
- ▶ 45° click feature
- ▶ Products can withstand upto 2 bar liquid pressure
- ▶ Hydrophobic filter prevents leakage while priming
- ▶ Rotating luer lock




REF	DESCRIPTION	SIZE - Tube Inner/outer, L(length)	PRIMING VOLUME(ml) Approx	Pack Count dispenser/carton
10548-01	3-way stopcock, Extension line, PVC-free	Ø 3/4.1 mm, L 10cm	0,089	25 / 100
10549-01	3-way stopcock, Extension line, PVC-free	Ø 3/4.1 mm, L 30cm	0,070	25 / 100
10550-01	3-way stopcock, Extension line, PVC-free	Ø 3/4.1 mm, L 50cm	0,066	25 / 100
10551-01	3-way stopcock, Extension line, PVC-free	Ø 3/4.1 mm, L 60cm	0,065	25 / 100
10552-01	3-way stopcock, Extension line, PVC-free	Ø 3/4.1 mm, L 90cm	0,066	25 / 100
10553-01	3-way stopcock, Extension line, PVC-free	Ø 3/4.1 mm, L 100cm	0,066	25 / 100
10713	3-way stopcock, Extension line, PVC-free	Ø 3/4.1 mm, L 15cm	0,086	25 / 100

Sales Unit: Transport Carton

## Regulations

The product complies with legal requirements for medical devices. EO sterilized according to valid version of ISO 11135 standard.



**REACH regulation:** The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).

**Biological evaluation:** EN ISO 10993-01

**Compliance with legal requirements:** Medical Device Regulation (EU) 2017/745, MDR class IIa

**Compliance with product standards:** EN 8536-9, EN ISO 80369-7

**Quality standards:** EN ISO 13485

**Label information and packaging compliance:** EN ISO 15223-1, EN ISO 20417

## Material

Product Part	Raw Material & Characteristics
Male luer cap	HDPE, white
3-way stopcock	Transparent channel PC, lipid resistant, white handle HDPE
Tube	Polyolefin alloy, PVC-free, soft, clear, kink resistant
Rotating luer lock	ABS & PC, transparent
Purge filter cap	HDPE, light blue

## Warnings & Recommendations

### Warnings

The use of this product is restricted to a qualified doctor or nurse.  
 Read the instructions before use. The device should be used in accordance with the instructions.  
 The device is guaranteed sterile and pyrogen-free if the individual package has not been opened or damaged. Do not use if package is damaged or if protective caps are loose or missing.  
 Do not over-tighten a luer beyond the friction fit as this may damage both the luer and another connecting device.  
 Do not use any instruments for clamping the tubing or luer tightening.  
 Re-use may lead to infection or other illness/injury.  
 Do not leave open ports/hubs exposed. Replace sterile end caps, as necessary.

### Storage Recommendation

Recommended storage in clean, dry space between 15 - 25 °C. Products should be protected from direct sunlight, other intensive light sources and ozone.

### Disposal Recommendation

Dispose in accordance with local regulations.

### Shelf Life

5 years

## Packing Information

Package	Material	Size ( L x W x H ) mm
Single Pack	Medical grade paper, PP/PE film	10548-01: 140 x 80
		10549-01: 235 x 100
		10550-01: 235 x 100
		10551-01: 235 x 100
		10552-01: 235 x 100
		10553-01: 235 x 100
		10713: 235 x 100
Dispenser Pack	PP	305 x 460
Transport Carton	TCF quality, Carton	10548-01: 310 x 270 x 192
		10549-01: 400 x 295 x 287
		10550-01: 400 x 295 x 287
		10551-01: 400 x 295 x 287
		10552-01: 400 x 295 x 287
		10553-01: 400 x 295 x 287
10713: 400 x 295 x 287		
Packing Tape	PVC-free, Transparent	N/A

## Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
10548-01	06438129117318	N/A	06438129521542
10549-01	06438129121551	N/A	06438129521559
10550-01	06438129121544	N/A	06438129521955
10551-01	06438129121964	N/A	06438129521962
10552-01	06438129121971	N/A	06438129521979
10553-01	06438129121988	N/A	06438129521986
10713	06438129125115	N/A	06438129525113

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.