


## evercare® Hip set standard

REF: 1210-02

evercare® Hip set standard, sterile is a draping set containing materials for hip operations. The component in size 180 x 180 cm and 150 x 240 cm are designed with adhesive tape with cut peel-papers in shorter length to ease the draping and allow sequenced covering.

### Key Features

- ▶ Absorbent patch
- ▶ Elastic stockinette
- ▶ Tube holders
- ▶ Adhesive OP Tapes
- ▶ Bandage
- ▶ Main drape complies to EN 13795-1 High performance critical product area requirements



REF	DESCRIPTION	SIZE	Pack Count dispenser/carton
1210-02	Hip set standard, sterile	See Material table	3 / 6

Sales Unit: Transport Carton

## Regulations

The product complies with legal requirements for medical devices. EO sterilized according to valid version of ISO 11135 standard.



### REACH regulation:

The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).

### Biological evaluation:

EN ISO 10993-01  
EN ISO 10993-5  
EN ISO 10993-7  
EN ISO 10993-10

### Compliance with legal requirements:

Medical Device Regulation (EU) 2017/745, MDR class Is

### Compliance with product standards:

EN13795-1

### Quality standards:

EN ISO 13485

### Label information and packaging compliance:

EN ISO 15223-1, EN ISO 20417, EN ISO 11607-1

## Material

Product Part	Raw Material & Characteristics
1 Mayo cover	78 x 145 cm, reinforced, Laminate: blue PP nonwoven/blue-white PE-film/blue PE-film
4 Hand towels	30 x 34 cm, Pulp, natural white, thread reinforced
1 Adhesive drape medium	180 x 180 cm, Laminate: blue PP nonwoven/blue PE-film/white PP nonwoven, blue absorbent patch 25 x 60 cm, PP nonwoven, tube holders
1 Elastic stockinette	30 x 120 cm, SEBS based thermoplastic elastomer/PES lining
1 Elastic support bandage	12 cm x 6 m PET/cotton
2 Adhesive drape sheets	75 x 90 cm, Laminate: blue PP nonwoven/blue PE-film/white PP nonwoven, blue absorbent patch 25 x 60 cm, PP nonwoven
3 Adhesive Op-tapes	9 x 50 cm, Nonwoven polyester, blue, adhesive coated
1 Adhesive split sheet	230 x 260 cm, Laminate: blue PP nonwoven/blue PE-film/white PP nonwoven, absorbent patch blue 140 x 75 cm, PP nonwoven, split 20 x 100, Tube holders
1 Adhesive drape large	150 x 240 cm, Laminate: blue PP nonwoven/blue PE-film/white PP nonwoven, blue absorbent patch 25 x 60 cm, PP nonwoven, tube holders
1 Instrument table cover	150 x 190 cm, Laminate: blue PP nonwoven/blue PE-film

## Warnings & Recommendations

<b>Warnings</b>	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur. Do not use if package is damaged.
<b>Storage Recommendation</b>	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
<b>Disposal Recommendation</b>	Dispose in accordance with local regulations.
<b>Shelf Life</b>	5 years

## Packing Information

Package	Material	Size ( L x W x H ) mm
<b>Single Pack</b>	Medical grade paper, PE/PET film	420 x 603
<b>Dispenser Pack</b>	TCF or ECF quality	485 x 285 x 287
<b>Transport Carton</b>	TCF or ECF quality	500 x 310 x 590
<b>Packing Tape</b>	Transparent	N/A

## Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
<b>1210-02</b>	06438129116717	06438129316711	06438129516715

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.